

PerCP/Cy5.5® Anti-Ly-6A/E (Sca-1) antibody [D7] ab93538

画像数 1

製品の概要

製品名	PerCP/Cy5.5® Anti-Ly-6A/E (Sca-1) antibody [D7]
製品の詳細	PerCP/Cy5.5® Rat monoclonal [D7] to Ly-6A/E (Sca-1)
由来種	Rat
標識	PerCP/Cy5.5®. Ex: 482nm, Em: 690nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Mouse
免疫原	Tissue, cells or virus corresponding to Mouse Sca1/ Ly6A/E. IL2 dependent T cell line CTL-L (Mouse).
ポジティブ・コントロール	C57Bl/6 Mouse splenocytes
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.20 Preservative: 0.09% Sodium azide Constituents: 0.87% Sodium chloride, 0.1% Gelatin, PBS
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	D7
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab93538の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		

追加情報

Flow Cyt: Use 0.125 - 0.25 µg for 10^5 - 10^8 cells in a 100 µl volume.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

機能	T-cell activation.
組織特異性	Widely expressed.
配列類似性	Contains 1 UPAR/Ly6 domain.
翻訳後修飾	O-glycosylated. Not N-glycosylated. Not phosphorylated.
細胞内局在	Cell membrane.

画像



Flow Cytometry - PerCP/Cy5.5® Anti-Ly-6A/E (Sca-1) antibody [D7] (ab93538)

Flow cytometry staining of C57Bl/6 splenocytes with 0.25 µg of PerCP-Cy5.5 Rat IgG2a Isotype Control (open histogram) or 0.125 µg of ab93538 (colored histogram). Cells in the lymphocyte gate were used for analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors