

PerCP/Cy5.5® Mouse IgG1, Kappa Monoclonal [MOPC-21] - Isotype Control ab157226

製品の概要

製品名	PerCP/Cy5.5® Mouse IgG1, Kappa Monoclonal [MOPC-21] - Isotype Control
標識	PerCP/Cy5.5®. Ex: 482nm, Em: 690nm
特異性	After screening on a variety of resting, activated, live and fixed rat and human tissues it was observed that ab157226 does not react specifically with rat and Human samples, and hence all the background that could be observed when working with ab157226 would be a result of general nonspecific interactions between an mouse IgG1 molecule and the respective sample under the particular conditions.
アプリケーション	適用あり: Flow Cyt
免疫原	ab157226 was generated against an undefined antigen.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	<p>pH: 7.4</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituents: 0.2% BSA, 99% PBS</p>
精製度	Size exclusion
特記事項 (精製)	The purified antibody is conjugated with tandem dye PerCP-Cy5.5 under optimum conditions. The conjugate is purified by size-exclusion chromatography.
ポリ/モノ	モノクローナル
クローン名	MOPC-21
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab157226の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent concentration. Use the same concentration of this isotype control antibody as the recommended working concentration of the antigen-specific antibody.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors