


Anti-Perilipin 2 antibody ab52356

★★★★★ **3 Abreviews** **31 References** 画像数 3

製品の概要

製品名	Anti-Perilipin 2 antibody
製品の詳細	Rabbit polyclonal to Perilipin-2
由来種	Rabbit
アプリケーション	適用あり: ICC, WB, IHC-P
種交差性	交差種: Mouse, Human 交差が予測される動物種: Rat 
免疫原	Synthetic peptide corresponding to Mouse Perilipin-2 (C terminal).
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.40 Constituents: 30% Glycerol, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab52356の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC		Use at an assay dependent concentration.
WB	★★★★★ (2)	Use a concentration of 2 µg/ml. Detects a band of approximately 48 kDa (predicted molecular weight: 48 kDa).
IHC-P		1/200.

ターゲット情報

機能	May be involved in development and maintenance of adipose tissue.
組織特異性	Milk lipid globules.
配列類似性	Belongs to the perilipin family.
翻訳後修飾	Acylated; primarily with C14, C16 and C18 fatty acids.
細胞内局在	Membrane.


画像



Anti-Perilipin 2 antibody (ab52356) at 2 µg/ml + mouse liver lysate


Predicted band size: 48 kDa

Observed band size: 48 kDa



Immunohistochemical analysis of mouse liver tissue staining Perilipin-2 using ab52356.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Perilipin 2 antibody (ab52356)



Ab52356 staining Perilipin-2 in HeLa cells at 1/50 dilution. Hoechst 33258 was used as a nuclear counterstain.

Immunocytochemistry - Anti-Perilipin 2 antibody (ab52356)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors