

Anti-PKC alpha antibody [133] ab11723

KO 評価済

リコンビナント

★★★★☆ 1 Abreviews 18 References 画像数 9

製品の概要

製品名	Anti-PKC alpha antibody [133]
製品の詳細	Mouse monoclonal [133] to PKC alpha
由来種	Mouse
特異性	ab11723 recognises the alpha isoform of PKC.
アプリケーション	適用あり: WB, IHC-P 適用なし: ICC or ICC/IF
種交差性	交差種: Mouse, Rat, Human
免疫原	Recombinant fragment. This information is proprietary to Abcam and/or its suppliers. (Peptide available as ab55672)
ポジティブ・コントロール	WB: HAP1, K562, C6, HeLa, NIH/3T3 and HEK-293 cell lysate. Human kidney tissue lysate. Mouse and rat brain and spleen tissue lysate. Recombinant human PKC alpha protein (Active) (ab55672). IHC-P: Human, mouse and rat retina tissue.
特記事項	<p>This product has switched from a hybridoma to recombinant production method on 25th July 2021.</p> <p>This antibody clone is manufactured by Abcam. If you require a custom buffer formulation or conjugation for your experiments, please contact orders@abcam.com.</p> <p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none">- High batch-to-batch consistency and reproducibility- Improved sensitivity and specificity- Long-term security of supply- Animal-free production <p>For more information see here.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.20 Preservative: 0.01% Sodium azide

	Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	133
ミエローマ	Sp2/0
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab11723の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB	★★★★★ (1)	1/500 - 1/1000. Predicted molecular weight: 76 kDa. Can be blocked with Recombinant human PKC alpha protein (Active) (ab55672) .
IHC-P		1/1 - 1/500. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.

追加情報 Is unsuitable for ICC or ICC/IF.

ターゲット情報

機能	This is a calcium-activated, phospholipid-dependent, serine- and threonine-specific enzyme. May play a role in cell motility by phosphorylating CSPG4. PKC is activated by diacylglycerol which in turn phosphorylates a range of cellular proteins. PKC also serves as the receptor for phorbol esters, a class of tumor promoters.
配列類似性	Belongs to the protein kinase superfamily. AGC Ser/Thr protein kinase family. PKC subfamily. Contains 1 AGC-kinase C-terminal domain. Contains 1 C2 domain. Contains 2 phorbol-ester/DAG-type zinc fingers. Contains 1 protein kinase domain.
細胞内局在	Cytoplasm. Cell membrane. Nucleus.

画像



Western blot - Anti-PKC alpha antibody [133]
(ab11723)

This image was generated from the hybridoma version of the product.

Lanes 1, 5 and 9: Wild-type HAP1 cell lysate (20 µg)

Lanes 2, 6 and 10: PKC alpha knockout HAP1 cell lysate (20 µg)

Lanes 3, 7 and 11: K562 cell lysate (20 µg)


Lanes 4, 8 and 12: HEK293 cell lysate (20 µg)

Lanes 1, 2, 3 and 4: Green signal from target – ab11723 observed at 77 kDa

Lanes 5, 6, 7 and 8: Red signal from loading control – **ab8245** observed at 37 kDa

Lanes 9, 10, 11 and 12: Merged (red and green) signal


ab11723 was shown to specifically react with PKC alpha in wild-type HAP1 cells. No band was observed when PKC alpha knockout samples were examined. Wild-type and PKC alpha knockout samples were subjected to SDS-PAGE. ab11723 and **ab8245** (loading control to GAPDH) were diluted 1/500 and 1/2000 respectively and incubated overnight at 4°C. Blots were developed with Goat anti-Mouse IgG H&L (IRDye® 800CW) preadsorbed (**ab216772**) and Goat Anti-Rabbit IgG H&L (IRDye® 680RD) preadsorbed (**ab216777**) secondary antibodies at 1/10,000 dilution for 1 hour at room temperature before imaging.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PKC alpha antibody [133]
(ab11723)

Immunohistochemical analysis of Paraffin-embedded human retina tissue labelling PKC alpha with ab11723 at 1/500 dilution, followed by ready to use secondary antibody Goat Anti-Mouse IgG H&L (HRP polymer) (**ab214879**). Positive staining on human retina is observed. Counter stained with Hematoxylin. Heat mediated antigen retrieval was performed using **ab93684** (Tris/EDTA buffer, pH 9.0).

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is ready to use Goat Anti-Mouse IgG H&L (HRP polymer) **ab214879**.



Western blot - Anti-PKC alpha antibody [133]
(ab11723)

This image was generated from the hybridoma version of the product.

Lane 1: Wild-type HAP1 cell lysate (20 µg)


Lane 2: PKC alpha knockout HAP1 cell lysate (20 µg)

Lane 3: K562 cell lysate (20 µg)

Lane 4: HEK293 cell lysate (20 µg)

Lanes 1 - 4: Merged signal (red and green). Green - ab11723 observed at 77 kDa. Red - loading control, **ab8245**, observed at 37 kDa.

This western blot image is a comparison between ab11723 and a competitor's top cited rabbit polyclonal antibody.



Western blot - Anti-PKC alpha antibody [133]
(ab11723)

All lanes : Anti-PKC alpha antibody [133] (ab11723) at 1/1000 dilution

Lane 1 : Recombinant human PKC alpha protein (Active)
(**ab55672**)

Lane 2 : Recombinant human PKC beta 1 protein (**ab60840**)

Lane 3 : Recombinant human PKC beta 2 protein (**ab60841**)

Lane 4 : Recombinant human PKC eta protein (**ab60849**)

Lane 5 : Recombinant human PKC epsilon protein (**ab60847**)

Lane 6 : Recombinant human PKC theta/PRKCQ protein
(**ab56641**)

Lane 7 : Recombinant human PKC mu/PKD protein (**ab60873**)

Lane 8 : Recombinant human PKC zeta protein (**ab60848**)

Lane 9 : Recombinant human PKC iota protein (**ab60850**)

Lane 10 : Recombinant human PKC gamma protein (**ab60842**)


Secondary

All lanes : Peroxidase-Conjugated Goat anti-Mouse IgG (H+L) at 1/5000 dilution

Predicted band size: 76 kDa

5% NFDN/TBST

All positive controls 10 ng



Western blot - Anti-PKC alpha antibody [133] (ab11723)

All lanes : Anti-PKC alpha antibody [133] (ab11723) at 1/1000 dilution

Lane 1 : Recombinant human PKC alpha protein (Active) ([ab55672](#))

Lane 2 : Recombinant human PKC delta protein ([ab60844](#))


Secondary

All lanes : Peroxidase-Conjugated Goat anti-Mouse IgG (H+L) at 1/5000 dilution

Predicted band size: 76 kDa

Exposure time: 180 seconds


5% NFDN/TBST



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PKC alpha antibody [133] (ab11723)

Immunohistochemical analysis of paraffin-embedded mouse retina tissue labelling PKC alpha with ab11723 at 1/500 dilution, followed by ready to use secondary antibody Goat Anti-Mouse IgG H&L (HRP polymer) ([ab214879](#)). Positive staining on mouse retina is observed. Counter stained with Hematoxylin. Heat mediated antigen retrieval was performed using [ab93684](#) (Tris/EDTA buffer, pH 9.0).

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is ready to use Goat Anti-Mouse IgG H&L (HRP polymer) [ab214879](#).



All lanes : Anti-PKC alpha antibody [133] (ab11723) at 1/1000 dilution

Lane 1 : Human kidney tissue lysate 20 ug

Lane 2 : Mouse brain tissue lysate 20 ug

Lane 3 : Mouse spleen tissue lysate 20 ug

Lane 4 : Rat brain tissue lysate 20 ug

Lane 5 : Rat spleen tissue lysate 20 ug


Secondary

All lanes : Peroxidase-Conjugated Goat anti-Mouse IgG (H+L) at 1/5000 dilution

Predicted band size: 76 kDa

Exposure time: 180 seconds

5% NFDm/TBST



All lanes : Anti-PKC alpha antibody [133] (ab11723) at 1/1000 dilution

Lane 1 : C6 (rat glial tumor glial cell), whole cell lysate 20 ug

Lane 2 : HeLa (human cervix adenocarcinoma epithelial cell), whole cell lysate 20 ug

Lane 3 : NIH/3T3 (mouse embryonic fibroblast), whole cell lysate 20 ug

Secondary

All lanes : Peroxidase-Conjugated Goat anti-Mouse IgG (H+L) at 1/5000 dilution


Predicted band size: 76 kDa

Exposure time

Lane 1: 136 seconds

Lane 2 and 3: 15 seconds

Fresh lysates were used in lane 2 and 3.



Immunohistochemical analysis of paraffin-embedded rat retina tissue labelling PKC alpha with ab11723 at 1/500 dilution, followed by ready to use secondary antibody Goat Anti-Mouse IgG H&L (HRP polymer) (**ab214879**). Positive staining on rat retina is observed. Counter stained with Hematoxylin. Heat mediated antigen retrieval was performed using **ab93684** (Tris/EDTA buffer, pH 9.0)

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is ready to use Goat Anti-Mouse IgG H&L (HRP polymer) **ab214879**.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PKC alpha antibody [133] (ab11723)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors