

# Anti-PKC theta/PRKCQ antibody [M217] ab76197

**1 References** [画像数 2](#)

### 製品の概要

製品名	Anti-PKC theta/PRKCQ antibody [M217]
製品の詳細	Mouse monoclonal [M217] to PKC theta/PRKCQ
由来種	Mouse
特異性	ab76197 is not predicted to react with other PKC family members due to low sequence homology.
アプリケーション	<b>適用あり:</b> ICC/IF, WB
種交差性	<b>交差種:</b> Rat, Human
免疫原	Recombinant fragment. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: Jurkat whole cell lysate ( <a href="#">ab7899</a> ). ICC/IF: PC-12 cells.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 50% Glycerol, PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	M217
アイソタイプ	IgG2a

### アプリケーション

## The Abpromise guarantee

**Abpromise保証は、次のテスト済みアプリケーションにおけるab76197の使用に適用されます**


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/100.
WB		1/1000. Predicted molecular weight: 82 kDa.

## ターゲット情報

機能	This is a calcium-independent, phospholipid-dependent, serine- and threonine-specific enzyme. Essential for T-cell receptor (TCR)-mediated T-cell activation, but is dispensable during TCR-dependent thymocyte development. Links the TCR signaling complex to the activation of NF-kappa-B in mature T lymphocytes. Required for interleukin-2 (IL2) production. PKC is activated by diacylglycerol which in turn phosphorylates a range of cellular proteins. PKC also serves as the receptor for phorbol esters, a class of tumor promoters.
組織特異性	Skeletal muscle, megakaryoblastic cells and platelets.
配列類似性	Belongs to the protein kinase superfamily. AGC Ser/Thr protein kinase family. PKC subfamily. Contains 1 AGC-kinase C-terminal domain. Contains 1 C2 domain. Contains 2 phorbol-ester/DAG-type zinc fingers. Contains 1 protein kinase domain.
ドメイン	The C1 domain, containing the phorbol ester/DAG-type region 1 (C1A) and 2 (C1B), is the diacylglycerol sensor and the C2 domain is a non-calcium binding domain.
翻訳後修飾	Autophosphorylation at Thr-219 is required for targeting to the TCR and cellular function of PKC upon antigen receptor ligation.

## 画像



**Lane 1 :** Anti-PKC theta/PRKCQ antibody [M217] (ab76197) at 1/250 dilution


**Lane 2 :** Anti-PKC theta/PRKCQ antibody [M217] (ab76197) at 1/1000 dilution

**Lane 1 :** Jurkat cell lysate.

**Lane 2 :** Jurkat cell lysate

**Predicted band size:** 82 kDa

**Observed band size:** 79 kDa



NGF-differentiated PC-12 cells (rat adrenal gland pheochromocytoma cell line) stained for PKC theta/PRKCQ using ab76197 at 1/100 dilution in ICC/IF. Secondary antibody is a Cy3 conjugate.

Immunocytochemistry/ Immunofluorescence - Anti-PKC theta/PRKCQ antibody [M217] (ab76197)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors