

Anti-Pneumolysin antibody [PLY-4] ab71810

★★★★★ [1 Abreviews](#) [15 References](#) [画像数 1](#)

製品の概要

製品名	Anti-Pneumolysin antibody [PLY-4]
製品の詳細	Mouse monoclonal [PLY-4] to Pneumolysin
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Species independent
免疫原	Recombinant full length protein
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.50 Constituent: 0.06% Tris buffered saline
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	PLY-4
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab71810の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB	★★★★★ (1)	Use a concentration of 3 µg/ml.

ターゲット情報

細胞内局在

Plasma membrane

画像



Coomassie Blue-stained gel (lane 1) and Western blots of purified recombinant PLY, 3 µg per lane, with mouse monoclonal antibodies PLY-4 (lane 2), PLY-5 (lane 3) and PLY-7 (lane 4); a PLY-irrelevant mouse monoclonal antibody 1.4G8.66 served as control (lane 5).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors