


Anti-PNPLA3 antibody ab87875

画像数 1

製品の概要

製品名	Anti-PNPLA3 antibody
製品の詳細	Goat polyclonal to PNPLA3
由来種	Goat
アプリケーション	適用あり: IHC-P, WB, ICC
種交差性	交差種: Mouse 交差が予測される動物種: Rat 
免疫原	Synthetic peptide: CVRKARSRNI GTLHPFFNIN KCIRDGLQES LPD, corresponding to amino acids 73-104 of Mouse PNPLA3 Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Mouse skin sections
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.1% Sodium azide Constituents: 0.1% BSA, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab87875の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
IHC-P		1/250.
WB		Use at an assay dependent concentration.
ICC		Use at an assay dependent concentration.

ターゲット情報

関連性 PNPLA3 is a multifunctional enzyme which has both triacylglycerol lipase and acylglycerol O-acyltransferase activities. It contains one patatin domain.

細胞内局在 Membrane; Single pass type II membrane protein

画像



Immunohistochemistry on paraffin sections of mouse skin using ab87875 at 1/250 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PNPLA3 antibody (ab87875)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors