

# Anti-PPP4R2 antibody ab70631

**2 References** [画像数 3](#)

### 製品の概要

---

製品名	Anti-PPP4R2 antibody
製品の詳細	Rabbit polyclonal to PPP4R2
由来種	Rabbit
アプリケーション	<b>適用あり:</b> IHC-P, WB, IP
種交差性	<b>交差種:</b> Mouse, Human <b>交差が予測される動物種:</b> Horse, Cow, Pig, Chimpanzee, Rhesus monkey, Gorilla, Orangutan



**免疫原** Synthetic peptide corresponding to a region between residue 370 and the C-terminus (residue 417) of Human PPP4R2 (NP\_777567.1)

**特記事項**

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.1% BSA, Tris buffered saline
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

---

## The Abpromise guarantee


**Abpromise保証は、次のテスト済みアプリケーションにおけるab70631の使用に適用されます**

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/200.
WB		1/2000 - 1/10000. Detects a band of approximately 65 kDa (predicted molecular weight: 47 kDa).
IP		Use at 2-5 µg/mg of lysate.

## ターゲット情報

## 画像



**All lanes :** Anti-PPP4R2 antibody (ab70631) at 0.1 µg/ml

**Lane 1 :** HeLa whole cell lysate

**Lane 2 :** HEK293T whole cell lysate


**Lane 3 :** Jurkat whole cell lysate

**Lane 4 :** TCMK whole cell lysate

**Lane 5 :** NIH 3T3 whole cell lysate


Lysates/proteins at 50 µg per lane.

**Predicted band size:** 47 kDa



Immunohistochemical analysis of human ovarian cancer tissue, labelling PPP4R2 with ab70631 at 1/200 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PPP4R2 antibody (ab70631)



Detection of Human PPP4R2 by Western Blot of Immunoprecipitate. ab70631 at 1µg/ml staining PPP4R2 in HeLa whole cell lysate immunoprecipitated using ab70631 at 3µg/mg lysate (1 mg/IP; 20% of IP loaded/lane). Detection: chemiluminescence with exposure time of 10 seconds.

Immunoprecipitation - Anti-PPP4R2 antibody (ab70631)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors