

Anti-Pseudomonas antibody ab68538

★★★★★ **3 Abreviews** **18 References**

製品の概要

製品名	Anti-Pseudomonas antibody
製品の詳細	Rabbit polyclonal to Pseudomonas
由来種	Rabbit
特異性	This antibody reacts with Pseudomonas in bacterial samples. It was raised against and validated with strain Boston 41501.
アプリケーション	適用あり: ICC/IF, IHC-Fr
種交差性	交差種: Pseudomonas
免疫原	Tissue, cells or virus corresponding to Pseudomonas. Whole cells of Pseudomonas aeruginosa.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.2 Preservative: 0.1% Sodium azide Constituent: 0.0268% PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、 次のテスト済みアプリケーションにおけるab68538の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ICC/IF	★★★★★ (1)	Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration. PubMed: 25375398

ターゲット情報

関連性 Pseudomonas are Gram-negative rod bacteria commonly found in soil, ground water, plants and animals. The most common form is Pseudomonas aeruginosa which is an opportunistic pathogen of both humans and plants.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors