

### Anti-PTEN antibody [EPR9941-2] ab170941

KO 評価済 リコンビナント RabMAb

★★★★☆ 1 Abreviews 47 References 画像数 10

#### 製品の概要

|              |   |
|--------------|---|
| 製品名          | Anti-PTEN antibody [EPR9941-2]  |
| 製品の詳細        | Rabbit monoclonal [EPR9941-2] to PTEN   |
| 由来種          | Rabbit  |
| アプリケーション     | 適用あり: Flow Cyt (Intra), IHC-P, WB   |
| 種交差性         | 交差種: Mouse, Rat, Human  |
| 免疫原          | Recombinant full length protein. This information is proprietary to Abcam and/or its suppliers.   |
| ポジティブ・コントロール | WB: HeLa, MCF7 and 293T cell lysates. IHC-P: Human breast carcinoma and colon tissues.  |
| 特記事項         | <p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> <li>- High batch-to-batch consistency and reproducibility</li> <li>- Improved sensitivity and specificity</li> <li>- Long-term security of supply</li> <li>- Animal-free production</li> </ul> <p>For more information <a href="#">see here</a>.</p> <p>Our RabMAb<sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <a href="#">RabMAb<sup>®</sup> patents</a>.</p> |

#### 製品の特性

|        |   |
|--------|---|
| 製品の状態  | Liquid  |
| 保存方法   | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. |
| バッファー  | Preservative: 0.01% Sodium azide<br>Constituents: 40% Glycerol, 59% PBS, 0.05% BSA  |
| 精製度    | Protein A purified  |
| ポリ/モノ  | モノクローナル   |
| クローン名  | EPR9941-2   |
| アイソタイプ | IgG   |

**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab170941の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション         | Abreviews | 特記事項   |
|------------------|-----------|--|
| Flow Cyt (Intra) |           | 1/120.   |
| IHC-P            |           | 1/50 - 1/100. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol. |
| WB               |           | 1/1000 - 1/5000. Predicted molecular weight: 47 kDa.   |

## ターゲット情報

|       |  |
|-------|--|
| 機能    | <p>Tumor suppressor. Acts as a dual-specificity protein phosphatase, dephosphorylating tyrosine-, serine- and threonine-phosphorylated proteins. Also acts as a lipid phosphatase, removing the phosphate in the D3 position of the inositol ring from phosphatidylinositol 3,4,5-trisphosphate, phosphatidylinositol 3,4-diphosphate, phosphatidylinositol 3-phosphate and inositol 1,3,4,5-tetrakisphosphate with order of substrate preference in vitro <math>\text{PtdIns}(3,4,5)\text{P}_3 &gt; \text{PtdIns}(3,4)\text{P}_2 &gt; \text{PtdIns}3\text{P} &gt; \text{Ins}(1,3,4,5)\text{P}_4</math>. The lipid phosphatase activity is critical for its tumor suppressor function. Antagonizes the PI3K-AKT/PKB signaling pathway by dephosphorylating phosphoinositides and thereby modulating cell cycle progression and cell survival. The unphosphorylated form cooperates with AIP1 to suppress AKT1 activation. Dephosphorylates tyrosine-phosphorylated focal adhesion kinase and inhibits cell migration and integrin-mediated cell spreading and focal adhesion formation. Plays a role as a key modulator of the AKT-mTOR signaling pathway controlling the tempo of the process of newborn neurons integration during adult neurogenesis, including correct neuron positioning, dendritic development and synapse formation. May be a negative regulator of insulin signaling and glucose metabolism in adipose tissue. The nuclear monoubiquitinated form possesses greater apoptotic potential, whereas the cytoplasmic nonubiquitinated form induces less tumor suppressive ability. In motile cells, suppresses the formation of lateral pseudopods and thereby promotes cell polarization and directed movement.</p> <p>Isoform alpha: Functional kinase, like isoform 1 it antagonizes the PI3K-AKT/PKB signaling pathway. Plays a role in mitochondrial energetic metabolism by promoting COX activity and ATP production, via collaboration with isoform 1 in increasing protein levels of PINK1.</p> |
| 組織特異性 | Expressed at a relatively high level in all adult tissues, including heart, brain, placenta, lung, liver, muscle, kidney and pancreas.   |
| 関連疾患  | <p>Cowden syndrome 1</p> <p>Lhermitte-Duclos disease</p> <p>Bannayan-Riley-Ruvalcaba syndrome</p> <p>Squamous cell carcinoma of the head and neck</p> <p>Endometrial cancer</p> <p>PTEN mutations are found in a subset of patients with Proteus syndrome, a genetically heterogeneous condition. The molecular diagnosis of PTEN mutation positive cases classifies Proteus syndrome patients as part of the PTEN hamartoma syndrome spectrum. As such, patients surviving the early years of Proteus syndrome are likely at a greater risk of developing malignancies.</p>   |

Glioma 2  
 VACTERL association with hydrocephalus  
 Prostate cancer  
 Macrocephaly/autism syndrome  
 A microdeletion of chromosome 10q23 involving BMPR1A and PTEN is a cause of chromosome 10q23 deletion syndrome, which shows overlapping features of the following three disorders: Bannayan-Zonana syndrome, Cowden disease and juvenile polyposis syndrome.

#### 配列類似性

Contains 1 C2 tensin-type domain.  
 Contains 1 phosphatase tensin-type domain.

#### ドメイン

The C2 domain binds phospholipid membranes in vitro in a Ca(2+)-independent manner; this binding is important for its tumor suppressor function.


#### 翻訳後修飾

Constitutively phosphorylated by CK2 under normal conditions. Phosphorylated in vitro by MAST1, MAST2, MAST3 and STK11. Phosphorylation results in an inhibited activity towards PIP3. Phosphorylation can both inhibit or promote PDZ-binding. Phosphorylation at Tyr-336 by FRK/PTK5 protects this protein from ubiquitin-mediated degradation probably by inhibiting its binding to NEDD4. Phosphorylation by ROCK1 is essential for its stability and activity. Phosphorylation by PLK3 promotes its stability and prevents its degradation by the proteasome. Monoubiquitinated; monoubiquitination is increased in presence of retinoic acid. Deubiquitinated by USP7; leading to its nuclear exclusion. Monoubiquitination of one of either Lys-13 and Lys-289 amino acid is sufficient to modulate PTEN compartmentalization. Ubiquitinated by XIAP/BIRC4.

#### 細胞内局在

Secreted. May be secreted via a classical signal peptide and reenter into cells with the help of a poly-Arg motif and Cytoplasm. Nucleus. Nucleus, PML body. Monoubiquitinated form is nuclear. Nonubiquitinated form is cytoplasmic. Colocalized with PML and USP7 in PML nuclear bodies. XIAP/BIRC4 promotes its nuclear localization.

#### 画像



**All lanes :** Anti-PTEN antibody [EPR9941-2] (ab170941) at 0.05 µg/ml (purified)

**Lane 1 :** HeLa (Human cervix adenocarcinoma epithelial cell) whole cell lysates

**Lane 2 :** 293T (Human embryonic kidney epithelial cell) whole cell lysates

**Lane 3 :** Rat brain lysates

**Lane 4 :** Mouse brain lysates


Lysates/proteins at 20 µg per lane.

#### Secondary

**All lanes :** Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution


**Predicted band size:** 47 kDa

Blocking and diluting buffer: 5% NFDM/TBST.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PTEN antibody  
[EPR9941-2] (ab170941)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Rat cerebrum tissue sections labeling PTEN with Purified ab170941 at 1:50 dilution (5.3 µg/ml). Heat mediated antigen retrieval was performed using **ab93684** (Tris/EDTA buffer, pH 9.0) ImmunoHistoProbe one step HRP Polymer (ready to use) was used as the secondary antibody. Negative control: PBS instead of the primary antibody. Hematoxylin was used as a counterstain




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PTEN antibody  
[EPR9941-2] (ab170941)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Mouse cerebrum tissue sections labeling PTEN with Purified ab170941 at 1:50 dilution (5.3 µg/ml). Heat mediated antigen retrieval was performed using **ab93684** (Tris/EDTA buffer, pH 9.0) ImmunoHistoProbe one step HRP Polymer (ready to use) was used as the secondary antibody. Negative control: PBS instead of the primary antibody. Hematoxylin was used as a counterstain



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PTEN antibody [EPR9941-2] (ab170941)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Human pancreas tissue sections labeling PTEN with Purified ab170941 at 1:50 dilution (5.3 µg/ml). Heat mediated antigen retrieval was performed using **ab93684** (Tris/EDTA buffer, pH 9.0) ImmunoHistoProbe one step HRP Polymer (ready to use) was used as the secondary antibody. Negative control: PBS instead of the primary antibody. Hematoxylin was used as a counterstain.




Western blot - Anti-PTEN antibody [EPR9941-2] (ab170941)

**Lane 1:** Wild-type HAP1 cell lysate (20 µg)

**Lane 2:** PTEN knockout HAP1 cell lysate (20 µg)


**Lanes 1 and 2:** Merged signal (red and green). Green - ab170941 observed at 47 kDa. Red - loading control, **ab8245**, observed at 37 kDa.

Unpurified ab170941 was shown to specifically react with PTEN in wild-type HAP1 cells. No band was observed when PTEN knockout samples were used. Wild-type and PTEN knockout samples were subjected to SDS-PAGE, ab170941 and **ab8245** (loading control to GAPDH) were diluted 1/500 and 1/1000 respectively and incubated overnight at 4°C. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed (**ab216773**) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed (**ab216776**) secondary antibodies at 1/10,000 dilution for 1 hr at room temperature before imaging.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PTEN antibody [EPR9941-2] (ab170941)

Immunohistochemical analysis of paraffin-embedded Human breast carcinoma tissue labeling PTEN with unpurified ab170941 at 1/50 dilution. Heat mediated antigen retrieval was performed using citrate buffer pH 6 before commencing with IHC staining protocol.



Flow Cytometry (Intracellular) - Anti-PTEN antibody [EPR9941-2] (ab170941)

Intracellular Flow Cytometry analysis of HeLa (human cervix adenocarcinoma) cells labeling PTEN with purified ab170941 at 1/120 dilution (10ug/ml) (red). Cells were fixed with 80% methanol and permeabilised with 0.1% Tween-20. A Goat anti rabbit IgG (Alexa Fluor® 488) ([ab150077](#)) (1/2000 dilution) was used as the secondary antibody. Rabbit monoclonal IgG (Black) ([ab172730](#)) was used as the isotype control, cells without incubation with primary antibody and secondary antibody (Blue) were used as the unlabeled control.



Western blot - Anti-PTEN antibody [EPR9941-2] (ab170941)

**All lanes :** Anti-PTEN antibody [EPR9941-2] (ab170941) at 1/1000 dilution (Unpurified)

**Lane 1 :** HeLa cell lysates


**Lane 2 :** MCF7 cell lysates

**Lane 3 :** 293T cell lysates

Lysates/proteins at 10 µg per lane.

**Predicted band size:** 47 kDa





Immunohistochemical analysis of paraffin-embedded Human colon tissue labeling PTEN with unpurified ab170941 at 1/50 dilution. Heat mediated antigen retrieval was performed using citrate buffer pH 6 before commencing with IHC staining protocol.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PTEN antibody  
[EPR9941-2] (ab170941)

#### Why choose a recombinant antibody?



**Research with confidence**  
Consistent and reproducible results



**Long-term and scalable supply**  
Recombinant technology



**Success from the first experiment**  
Confirmed specificity



**Ethical standards compliant**  
Animal-free production

Anti-PTEN antibody [EPR9941-2] (ab170941)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors