

# Anti-PYGO2 antibody [PYGAD53A] ab51367

画像数 1

### 製品の概要

製品名	Anti-PYGO2 antibody [PYGAD53A]
製品の詳細	Mouse monoclonal [PYGAD53A] to PYGO2
由来種	Mouse
アプリケーション	<b>適用あり:</b> WB
免疫原	Recombinant fragment corresponding to Human PYGO2.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.05% Sodium azide</p> <p>Constituents: 1% BSA, PBS</p>
精製度	Protein G purified
特記事項 (精製)	Filtered through a 0.22 µm membrane.
ポリ/モノ	モノクローナル
クローン名	PYGAD53A
アイソタイプ	IgG1

### アプリケーション

The Abpromise guarantee

**Abpromise保証は、次のテスト済みアプリケーションにおけるab51367の使用に適用されます**


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration. Predicted molecular weight: 41 kDa.

## ターゲット情報

機能	Involved in signal transduction through the Wnt pathway.
配列類似性	Contains 1 PHD-type zinc finger.
細胞内局在	Nucleus.

## 画像



**Lane 1 :** Molecular weight markers

**Lane 2 :** Anti-PYGO2 antibody [PYGAD53A] (ab51367)

**Lane 2 :** immunizing peptide

**Predicted band size:** 41 kDa

**Observed band size:** 35,38 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors