

### Anti-RAB43 antibody [5G4] ab58030

[4 References](#) [画像数 2](#)

#### 製品の概要

製品名	Anti-RAB43 antibody [5G4]
製品の詳細	Mouse monoclonal [5G4] to RAB43
由来種	Mouse
アプリケーション	<b>適用あり:</b> WB, IHC-P
種交差性	<b>交差種:</b> Human
免疫原	Recombinant fragment corresponding to Human RAB43 aa 113-211. Sequence: IEDVRKYAGSNIVQLLIGNKSDLSELREVSLAEQSLAEHYD ILCAIETS AKDSSNVEEAFLRVATELIMRHGGPLFSEKSPDHIQLNSKDI GEGWGC

 [Run BLAST with](#)

 [Run BLAST with](#)

#### 特記事項

This product was changed from ascites to tissue culture supernatant on 11<sup>th</sup> June 2019. Please note that the dilutions may need to be adjusted accordingly. If you have any questions, please do not hesitate to contact our scientific support team.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.4
精製度	Tissue culture supernatant
ポリ/モノ	モノクローナル

クローン名	5G4
アイソタイプ	IgG2a
軽鎖の種類	kappa

## アプリケーション


**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab58030の使用に適用されず  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration. Predicted molecular weight: 23 kDa.
IHC-P		Use at an assay dependent concentration.

## ターゲット情報

組織特異性	Widely expressed in brain, testis, lung, heart, ovary, colon, kidney, uterus and spleen but not in liver.
配列類似性	Belongs to the small GTPase superfamily. Rab family.
細胞内局在	Cell membrane.


## 画像



RAB43 antibody (ab58030) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human salivary gland.

This image was generated using the ascites version of the product.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-RAB43 antibody [5G4] (ab58030)



RAB43 antibody (ab58030) at 1ug/lane + HepG2 cell lysate at 25ug/lane.

This image was generated using the ascites version of the product.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors