

Anti-RAP1/RABGEF1 antibody ab113480

3 References [画像数 1](#)

製品の概要

製品名	Anti-RAP1/RABGEF1 antibody
製品の詳細	Rabbit polyclonal to RAP1/RABGEF1
由来種	Rabbit
特異性	At least four isoforms of RAP1 are known to exist; ab113480 will recognize all four. ab113480 is predicted to not cross-react with other LRRFIP family members.
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to Human RAP1/RABGEF1 aa 1-100 (internal sequence).
ポジティブ・コントロール	Rat brain tissue lysate
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C long term.
バッファー	<p>pH: 7.2</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: 99.98% PBS</p>
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab113480の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.5 - 1 µg/ml. Predicted molecular weight: 79 kDa.

ターゲット情報

機能	Rab effector protein acting as linker between gamma-adaptin, RAB4A or RAB5A. Involved in endocytic membrane fusion and membrane trafficking of recycling endosomes. Stimulates nucleotide exchange on RAB5A. Can act as a ubiquitin ligase.
配列類似性	Contains 1 A20-type zinc finger. Contains 1 VPS9 domain.
翻訳後修飾	Monoubiquitinated.
細胞内局在	Cytoplasm. Early endosome. Recycling endosome.

画像



Lane 1 : Anti-RAP1/RABGEF1 antibody (ab113480) at 1 µg/ml

Lane 2 : Anti-RAP1/RABGEF1 antibody (ab113480) at 2 µg/ml

All lanes : Human liver tissue lysate

Predicted band size: 79 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors