

Anti-RASSF10 antibody ab113105

4 References [画像数 2](#)

製品の概要

製品名	Anti-RASSF10 antibody
製品の詳細	Rabbit polyclonal to RASSF10
由来種	Rabbit
特異性	ab113105 will not cross-react with other RAS association domain-containing family members.
アプリケーション	適用あり: IHC-P, WB
種交差性	交差種: Human
免疫原	Synthetic peptide containing a 17 amino acid sequence from the N terminal region of Human RASSF10 (NP_001073990).
ポジティブ・コントロール	Human lung tissue lysate.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	<p>pH: 7.2</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: 99.98% PBS</p>
精製度	Immunogen affinity purified
特記事項(精製)	ab113105 is purified by affinity chromatography via a peptide column.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab113105の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 5 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 57 kDa.

ターゲット情報


配列類似性 Contains 1 Ras-associating domain.

画像



Anti-RASSF10 antibody (ab113105) at 1 µg/ml + Human lung tissue lysate

Predicted band size: 57 kDa



Immunohistochemistry of RASSF10 in human lung tissue with ab113105 at 5 ug/mL.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors