

# Anti-Red Blood Cells antibody ab34858

## 6 References

### 製品の概要

製品名	Anti-Red Blood Cells antibody
製品の詳細	Rabbit polyclonal to Red Blood Cells
由来種	Rabbit
アプリケーション	<b>適用あり:</b> Agglutination
種交差性	<b>交差種:</b> Human
免疫原	Tissue/ cell preparation: Human washed pooled Red Blood Cells.
特記事項	<p>ab34858 is whole antiserum. For your information, we provide you the total protein concentration which is 90 mg/ml. Since the antibody is unpurified, we can only estimate its concentration at 1–10 mg/ml.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Sodium azide Constituents: 0.88% Sodium chloride, 0.42% Potassium phosphate
特記事項 (精製)	Ab34858 was prepared from polyspecific antiserum by a delipidation and defibrination.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab34858の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Agglutination		Use at an assay dependent concentration. An agglutination titer against serial dilutions of washed pooled cells on a titer plate of 1:64 was observed for this product.

## ターゲット情報

**関連性**      Red blood cells are the most common type of blood cell and are the vertebrate body's principal means of delivering oxygen from the lungs or gills to body tissues via the blood.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors