

Anti-Reovirus p17 protein antibody ab26810

画像数 1

製品の概要

製品名	Anti-Reovirus p17 protein antibody
製品の詳細	Rabbit polyclonal to Reovirus p17 protein
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Reovirus sp.
免疫原	Synthetic peptide corresponding to Reovirus p17 protein. Database link: Q99213

特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab26810の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/500. Predicted molecular weight: 17 kDa.
ELISA		1/50000.

ターゲット情報

関連性 The biological role of Reovirus p17 protein is fully unknown so far. It is possible that Reovirus p17 protein possesses a growth retardation function.

画像



Western blot analysis of Reovirus p17 protein.

All lanes : Anti-Reovirus p17 protein antibody (ab26810) at 1/500 dilution

Lane 1: p17 recombinant protein

Lane 2: pcDNA3.1

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors