

Anti-RPP30 antibody [2931D5a] ab53744

画像数 2

製品の概要

製品名	Anti-RPP30 antibody [2931D5a]
製品の詳細	Mouse monoclonal [2931D5a] to RPP30
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Human
免疫原	Recombinant fragment (Human)
ポジティブ・コントロール	HeLa whole cell lysate NIH3T3 whole cell lysate
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.05% Sodium azide</p> <p>Constituents: 1% BSA, PBS</p>
精製度	Protein G purified
特記事項 (精製)	Purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG depleted (approximately 95%) fetal bovine serum.
ポリ/モノ	モノクローナル
クローン名	2931D5a
アイソタイプ	IgG1


アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab53744の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/100. Detects a band of approximately 30 kDa (predicted molecular weight: 30 kDa).

ターゲット情報

画像



Western blot - Anti-RPP30 antibody [2931D5a]
(ab53744)

Anti-RPP30 antibody [2931D5a] (ab53744) at 1/100 dilution +
HeLa whole cell lysate at 25 µg


Secondary

anti Mouse IgG at 1/2500 dilution

Developed using the ECL technique.

Predicted band size: 30 kDa

Observed band size: 30 kDa



Western blot - Anti-RPP30 antibody [2931D5a]
(ab53744)

Anti-RPP30 antibody [2931D5a] (ab53744) at 1/100 dilution +
NIH3T3 whole cell lysate at 25 µg

Secondary

anti Mouse IgG at 1/2500 dilution

Developed using the ECL technique.

Predicted band size: 30 kDa

Observed band size: 30 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors