

# Anti-S100A9 + Calprotectin (S100A8/A9 complex) antibody [MAC387] - Low endotoxin, Azide free ab80084

**4 References** [画像数 1](#)

### 製品の概要

製品名	Anti-S100A9 + Calprotectin (S100A8/A9 complex) antibody [MAC387] - Low endotoxin, Azide free
製品の詳細	Mouse monoclonal [MAC387] to S100A9 + Calprotectin (S100A8/A9 complex) - Low endotoxin, Azide free
由来種	Mouse
アプリケーション	<b>適用あり:</b> Flow Cyt (Intra)
種交差性	<b>交差種:</b> Human
免疫原	Tissue, cells or virus corresponding to Human S100A9 + Calprotectin (S100A8/A9 complex). Solubilized human monocyte
ポジティブ・コントロール	Flow Cyt (Intra): Human peripheral blood granulocytes.
特記事項	<p>Endotoxin Level &lt;0.01EU/ug.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性


製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
キャリア・フリー	はい
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	MAC387
ミエローマ	NS1
アイソタイプ	IgG1

## アプリケーション

**The Abpromise guarantee** **Abpromise保証は、次のテスト済みアプリケーションにおけるab80084の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Flow Cyt (Intra)		1/50 - 1/200. Membrane permeabilisation is required for this application. <b>ab170190</b> - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

## 画像



**ab80085** at 1/50-1/200 dilution staining S100A9 + Calprotectin in human peripheral blood granulocytes by intracellular flow cytometry.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors