





Anti-SAMD12 antibody ab121831

画像数 6

製品の概要

製品名	Anti-SAMD12 antibody
製品の詳細	Rabbit polyclonal to SAMD12
由来種	Rabbit
アプリケーション	適用あり: ICC/IF, IHC-P, WB
種交差性	交差種: Human
免疫原	antigen sequence: LHCGLNPRGI DHPAHAEGIK LQIEGEGVES QSIKNKNFQK VPDQKGTPKR LQAEAETAKS ATVKLSKPVA LWTQQDVCKW LKKHCPNQYQ IYSESF, corresponding to amino acids 6-101 of Human SAMD12 (UniProt ID:Q8N8I0). <div>  Run BLAST with   Run BLAST with  </div>
ポジティブ・コントロール	Human brain tissue
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab121831の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 0.25 - 2 µg/ml.
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		Use a concentration of 0.04 - 0.4 µg/ml.

ターゲット情報


配列類似性 Contains 1 SAM (sterile alpha motif) domain.

画像




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SAMD12 antibody (ab121831)

Immunohistochemical analysis of human cerebral cortex tissue labelling SAMD12 using ab121831 at a 1/200 dilution. Cytoplasmic staining observed in human neuronal cells.




Immunohistochemical analysis of human pancreas tissue labelling SAMD12 using ab121831 at a 1/200 dilution. Cytoplasmic staining observed in human islets of Langerhans.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SAMD12 antibody (ab121831)




Immunohistochemical analysis of human fallopian tube tissue labelling SAMD12 using ab121831 at a 1/200 dilution. Cytoplasmic staining observed in human glandular cells.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SAMD12 antibody (ab121831)




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SAMD12 antibody (ab121831)

Immunohistochemical analysis of human colon tissue labelling SAMD12 using ab121831 at a 1/200 dilution. Cytoplasmic staining observed in human glandular cells.



Immunocytochemistry/ Immunofluorescence - Anti-SAMD12 antibody (ab121831)

Immunofluorescent staining of Human cell line U-251MG shows positivity in cytoplasm and golgi apparatus. Recommended concentration of ab121831 1-4 µg/ml. Cells treated with PFA/Triton X-100.



Western blot - Anti-SAMD12 antibody (ab121831)

All lanes : Anti-SAMD12 antibody (ab121831) at 1/250 dilution

Lane 1 : Negative control (vector only transfected HEK293T lysate)

Lane 2 : Over expression Lysate (co-expressed with a C terminal myc-DDK tag (~3.1 kDa) in mammalian HEK293T cells)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors