

### Anti-SCGB3A1 antibody ab92350

★★★★★ [1 Abreviews](#)

#### 製品の概要

製品名 Anti-SCGB3A1 antibody

製品の詳細 Goat polyclonal to SCGB3A1

① This product is a **fast track antibody**. It has been affinity purified and shows high titre values against the immunizing peptide by ELISA.

[Read the terms of use »](#)

由来種 Goat

種交差性

交差が予測される動物種: Human 

免疫原 Synthetic peptide:

NHLIEGSQKCVAE

, corresponding to internal sequence amino acids 69 - 81 of Human SCGB3A1 (NP\_443095.2)

 [Run BLAST with](#)

 [Run BLAST with](#)

特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

#### 製品の特性

製品の状態 Liquid

保存方法 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

バッファー

pH: 7.30

Preservative: 0.02% Sodium azide

Constituents: Tris buffered saline, 0.5% BSA

精製度

Immunogen affinity purified

特記事項 (精製)

ab92350 is purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

ポリ/モノ  
アイソタイプ

ポリクローナル  
IgG

#### アプリケーション

Fast track antibodies constitute a diverse group of products that have been released to accelerate your research, but are not yet fully characterized. They have all been affinity purified and show high titre values against the immunizing peptide (by ELISA).

#### **Fast track terms of use**

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration. Predicted molecular weight: 10 kDa.  A 1 hour primary incubation is recommended for this product. Preliminary experiments gave an approx 33kDa band in Human Breast cancer and normal Lung lysates after 0.5µg/ml antibody staining. Please note that currently we cannot find an explanation for this.

#### ターゲット情報

**機能** Potential growth inhibitory cytokine.

**組織特異性** Highly expressed in breast tissues. Absent in breast cancer cell lines.

**配列類似性** Belongs to the secretoglobin family. UGRP subfamily.

**細胞内局在** Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors