



### Anti-Shiga Toxin A subunit antibody [1C6A6G9] ab101839

★★★★★ **1 Abreviews**

#### 製品の概要

製品名	Anti-Shiga Toxin A subunit antibody [1C6A6G9]
製品の詳細	Mouse monoclonal [1C6A6G9] to Shiga Toxin A subunit
由来種	Mouse
アプリケーション	<b>適用あり:</b> ELISA, WB
種交差性	<b>交差種:</b> Escherichia coli
免疫原	Synthetic peptide (REFRQALSETAPVYT ) corresponding to amino acids 198-212 of Shiga Toxin A subunit.  <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a>

#### 特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Affinity purified
ポリ/モノ	モノクローナル
クローン名	1C6A6G9
アイソタイプ	IgG1

#### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、**次のテスト済みアプリケーションにおけるab101839の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
WB	★★★★☆ (1)	Use at an assay dependent concentration.

## ターゲット情報

**関連性**      Shiga toxins are a family of related toxins with two major groups, Stx1 and Stx2, whose genes are considered to be part of the genome of lambdoid prophages. The most common sources for Shiga toxin are the bacteria *S. dysenteriae* and the Shigatoxigenic group of *Escherichia coli* (STEC), which includes serotype O157:H7 and other enterohemorrhagic *E. coli* (EHEC).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors