

# Anti-SLITRK4 antibody ab67315

**1 References**   [画像数 1](#)

### 製品の概要

製品名	Anti-SLITRK4 antibody
製品の詳細	Rabbit polyclonal to SLITRK4
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Mouse
免疫原	A 13 amino acid peptide from near the C-terminus of Human SLITRK4.
ポジティブ・コントロール	WB: mouse brain tissue lysate
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション


**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab67315の使用に適用されます**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.5 - 1 µg/ml. Detects a band of approximately 94 kDa (predicted molecular weight: 94 kDa).

## ターゲット情報

機能	Suppresses neurite outgrowth.
組織特異性	Expressed in the cerebral cortex of the brain and at higher levels in some astrocytic brain tumors such as astrocytomas, glioblastomas and primitive neuroectodermal tumors.
配列類似性	Belongs to the SLITRK family. Contains 12 LRR (leucine-rich) repeats. Contains 2 LRRCT domains. Contains 1 LRRNT domain.
細胞内局在	Membrane.

## 画像



**Lane 1 :** Anti-SLITRK4 antibody (ab67315) at 0.5 µg/ml

**Lane 2 :** Anti-SLITRK4 antibody (ab67315) at 1 µg/ml

**All lanes :** mouse brain tissue lysate

Lysates/proteins at 15 µg per lane.

**Predicted band size:** 94 kDa

**Observed band size:** 94 kDa

**Additional bands at:** 120 kDa. We are unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors