


Anti-SOX2 antibody [20G5] ab171380

★★★★★ [1 Abreviews](#) [48 References](#) [画像数 8](#)

製品の概要

製品名	Anti-SOX2 antibody [20G5]
製品の詳細	Mouse monoclonal [20G5] to SOX2
由来種	Mouse
アプリケーション	適用あり: IHC-P, IP, WB, Flow Cyt, ICC/IF
種交差性	交差種: Mouse, Human 交差が予測される動物種: Sheep, Chicken, Xenopus laevis, Zebrafish, Xenopus tropicalis 
免疫原	Recombinant full length protein corresponding to Human SOX2 aa 1-317. Expressed in bacteria. Database link: 6657
ポジティブ・コントロール	IHC-P: Human lung squamous carcinoma tissue. Mouse esophagus tissue. ICC/IF: H9 embryonic stem cells. HEL 11.4 induced IPS cells. IP: NCCIT whole cell lysate. WB: NCCIT and NTERRA cell lysate. Flow Cytometry: H9 embryonic stem cells. HEL 11.4 induced IPS cells.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.05% Sodium azide Constituents: PBS, 30% Glycerol, 0.1% BSA
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	20G5
アイソタイプ	IgG1

アプリケーション


The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab171380の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P	★★★★★ (1)	1/20 - 1/200.
IP		Use at 5 µg/mg of lysate.
WB		1/1000 - 1/2000. Predicted molecular weight: 34 kDa.
Flow Cyt		1/100. ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.
ICC/IF		1/200.

ターゲット情報

機能	Transcription factor that forms a trimeric complex with OCT4 on DNA and controls the expression of a number of genes involved in embryonic development such as YES1, FGF4, UTF1 and ZFP206 (By similarity). Critical for early embryogenesis and for embryonic stem cell pluripotency.
関連疾患	Defects in SOX2 are the cause of microphthalmia syndromic type 3 (MCOPS3) [MIM:206900]. Microphthalmia is a clinically heterogeneous disorder of eye formation, ranging from small size of a single eye to complete bilateral absence of ocular tissues (anophthalmia). In many cases, microphthalmia/anophthalmia occurs in association with syndromes that include non-ocular abnormalities. MCOPS3 is characterized by the rare association of malformations including uni- or bilateral anophthalmia or microphthalmia, and esophageal atresia with trachoesophageal fistula.
配列類似性	Contains 1 HMG box DNA-binding domain.
翻訳後修飾	Sumoylation inhibits binding on DNA and negatively regulates the FGF4 transactivation.
細胞内局在	Nucleus.

画像



All lanes : Anti-SOX2 antibody [20G5] (ab171380) at 1/500 dilution

Lane 1 : NCCIT (Human pluripotent embryonic carcinoma cell line) cell lysate

Lane 2 : NTERRA cell lysate

Lane 3 : HeLa (Human epithelial cell line from cervix adenocarcinoma) cell lysate


Lysates/proteins at 25 µg per lane.

Secondary

All lanes : mouse IgG-HRP at 1/10000 dilution

Developed using the ECL technique.

Predicted band size: 34 kDa




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SOX2 antibody [20G5] (ab171380)

Image from Lu L. et al., Sci Rep. 2018 Aug 14;8(1):12083. Fig2c. doi: 10.1038/s41598-018-29441-3. Reproduced under the Creative Commons license <http://creativecommons.org/licenses/by/4.0/>.

REST expression in *N-hREST* mouse brains correlates with stemness in embryonic neural stem cells. Immunofluorescence analysis of E18.5 *N-hREST* and *LSL-hREST* control littermate mouse brains with antibodies against REST (using an antibody that preferentially recognizes hREST over mouse REST) and SOX2 (using ab171380).


Mice were anesthetized and perfused with phosphate-buffered saline followed by 4% paraformaldehyde (PFA). Brain tissues were then dissected and fixed in 4% PFA overnight at 4 °C. Fixed brain tissues were processed for paraffin embedding and then cut into 5-µm sections.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SOX2 antibody [20G5] (ab171380)

Immunohistochemistry analysis of SOX2 showing staining in the nucleus of paraffin-treated human lung squamous carcinoma (right) compared with a negative control without primary antibody (left).


To expose target proteins, antigen retrieval was performed using 10 mM sodium citrate (pH 6.0), microwaved for 8-15 minutes. Following antigen retrieval, tissues were blocked in 3% H₂O₂-methanol for 15 minutes at room temperature, washed with ddH₂O and PBS, and then probed with a SOX2 monoclonal antibody (ab171380) diluted by 3% BSA-PBS at a dilution of 1:200 overnight at 4°C in a humidified chamber. Tissues were washed extensively in PBST and detection was performed using an HRP-conjugated secondary antibody followed by colorimetric detection using a DAB kit. Tissues were counterstained with hematoxylin and dehydrated with ethanol and xylene to prep for mounting.



Immunocytochemistry/ Immunofluorescence - Anti-SOX2 antibody [20G5] (ab171380)


Immunofluorescence analysis of formaldehyde-fixed H9 embryonic stem cells, labeling SOX2 using ab171380 (left panel) at a 1/200 dilution overnight.

DAPI was used to stain the cell nuclei (central panel). Slides were washed with PBS and incubated with a fluorescein-conjugated secondary antibody at a 1/100 dilution.




Flow Cytometry - Anti-SOX2 antibody [20G5] (ab171380)

Flow cytometry analysis of H9 embryonic stem cells labeling SOX2 (blue histogram), using ab171380 at a 1/100 dilution, or a mouse IgG (green histogram) at a 1/100 dilution. A fluorescein-conjugated secondary antibody at a 1/200 dilution was used for the analysis.



Immunocytochemistry/ Immunofluorescence - Anti-SOX2 antibody [20G5] (ab171380)


Immunofluorescence analysis of formaldehyde-fixed HEL 11.4 induced IPS cells, labeling SOX2 using ab171380 (left panel) at a 1/200 dilution overnight. DAPI was used to stain the cell nuclei (central panel). Slides were washed with PBS and incubated with a fluorescein-conjugated secondary antibody at a 1/100 dilution.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SOX2 antibody [20G5] (ab171380)

Immunohistochemistry analysis of SOX2 showing staining in the nucleus of paraffin-treated mouse esophagus tissue (right) compared with a negative control without primary antibody (left).

To expose target proteins, antigen retrieval was performed using 10 mM sodium citrate (pH 6.0), microwaved for 8-15 minutes. Following antigen retrieval, tissues were blocked in 3% H₂O₂-methanol for 15 minutes at room temperature, washed with ddH₂O and PBS, and then probed with a SOX2 monoclonal antibody (ab171380) diluted by 3% BSA-PBS at a dilution of 1:20 overnight at 4°C in a humidified chamber. Tissues were washed extensively in PBST and detection was performed using an HRP-conjugated secondary antibody followed by colorimetric detection using a DAB kit. Tissues were counterstained with hematoxylin and dehydrated with ethanol and xylene to prep for mounting.



Flow Cytometry - Anti-SOX2 antibody [20G5] (ab171380)

Flow cytometry analysis of HEL 11.4 induced IPS cells labeling SOX2 (blue histogram), using ab171380 at a 1/100 dilution, or a mouse IgG (green histogram) at a 1/100 dilution. A fluorescein-conjugated secondary antibody at a 1/200 dilution was used for the analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors