

Anti-SPNS2 antibody ab82629

2 References [画像数 1](#)

製品の概要

製品名	Anti-SPNS2 antibody
製品の詳細	Rabbit polyclonal to SPNS2
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Mouse, Rat, Human, Drosophila melanogaster, Monkey
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	IHC: rat ovary tissue
特記事項	Storage of very dilute antibody solutions is not recommended The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium azide Stabilization buffer: contains Phosphate buffer; stabilizing proteins; cryogenic protectant and antibody stabilizer.
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab82629の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/150.

ターゲット情報

機能	Sphingolipid transporter required for migration of myocardial precursors. Transports sphingosine 1-phosphate (S1P), a secreted lipid mediator that plays critical roles in cardiovascular, immunological, and neural development and function. Mediates the export of S1P from cells in the extraembryonic yolk syncytial layer (YSL), thereby regulating myocardial precursor migration.
配列類似性	Belongs to the major facilitator superfamily. Spinster (TC 2.A.1.49) family.
細胞内局在	Membrane.

画像



Immunohistochemistry analysis of formalin-fixed paraffin-embedded rat ovary tissue staining SPNS2 with ab82629 at 1/50 dilution with DAB and Hematoxylin QS counterstaining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SPNS2 antibody (ab82629)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors