

Anti-ST2 antibody - BSA free ab25877

★★★★★ [2 Abreviews](#) [25 References](#) [画像数 6](#)

製品の概要

製品名	Anti-ST2 antibody - BSA free
製品の詳細	Rabbit polyclonal to ST2 - BSA free
由来種	Rabbit
特異性	ab25877 reacts all three known ST2 isoforms.
アプリケーション	適用あり: IHC-P, ICC/IF, WB
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide corresponding to Mouse ST2 aa 40-90 (N terminal). (Peptide available as ab39855)
ポジティブ・コントロール	IHC-P: Human lung cancer, Human liver and Mouse kidney tissues. WB: Hela, HepG2, 3T3/NIH, YB2/0 whole cell lysates. ICC/IF: Mouse kidney cells.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C. Avoid freeze / thaw cycle. Stable for 12 months at -20°C.
バッファー	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
特記事項(精製)	ST2 Antibody is affinity chromatography purified via peptide column.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab25877の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P	★★★★★ (2)	Use a concentration of 2 µg/ml.
ICC/IF		Use a concentration of 20 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 40 kDa. ab25877 at 2ug/ml. 1 hour incubation at room temperature in 5% NFD/MTBST. Secondary: Goat Anti-Rabbit IgG HRP conjugate at 1:10000 dilution.

ターゲット情報


機能	Receptor for interleukin-33 (IL-33), its stimulation recruits MYD88, IRAK1, IRAK4, and TRAF6, followed by phosphorylation of MAPK3/ERK1 and/or MAPK1/ERK2, MAPK14, and MAPK8. Possibly involved in helper T-cell function.
組織特異性	Highly expressed in kidney, lung, placenta, stomach, skeletal muscle, colon and small intestine. Isoform A is prevalently expressed in the lung, testis, placenta, stomach and colon. Isoform B is more abundant in the brain, kidney and the liver. Isoform C is not detected in brain, heart, liver, kidney and skeletal muscle.
配列類似性	Belongs to the interleukin-1 receptor family. Contains 3 Ig-like C2-type (immunoglobulin-like) domains. Contains 1 TIR domain.
細胞内局在	Secreted and Cell membrane.

画像




Human lung cancer tissue stained for ST2 using ab25877 at 5 µg/ml in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ST2 antibody - BSA free (ab25877)




Immunohistochemistry (Formalin-fixed paraffin embedded sections) of human lung cancer tissue labeling ST2 with Anti-ST2 antibody (ab25877) at 5µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ST2 antibody - BSA free (ab25877)




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ST2 antibody - BSA free (ab25877)

Immunohistochemistry of ST2 in human liver tissue with ST2 antibody at 5 µg/ml.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ST2 antibody - BSA free (ab25877)

ab25877 at 2µg/ml staining ST2 in mouse kidney tissue by IHC



Western blot - Anti-ST2 antibody - BSA free (ab25877)

All lanes : Anti-ST2 antibody - BSA free (ab25877) at 0.5 µg/ml

Lane 1 : HeLa (Human epithelial cell line from cervix adenocarcinoma) whole cell lysate

Lane 2 : HepG2 (Human liver hepatocellular carcinoma cell line) whole cell lysate

Lane 3 : NIH/3T3 (Mouse embryonic fibroblast cell line) whole cell lysate

Lane 4 : YB2/O whole cell lysate

Lysates/proteins at 15 µg per lane.

Secondary

All lanes : Rabbit IgG antibody (HRP) at 1/10000 dilution

Predicted band size: 40 kDa

Observed band size: 70 kDa

Exposure time: 1 minute

Western blot analysis labeling ST2 with ab25877 at 0.5 µg/mL in 5% NFDM/TBST, followed by secondary antibody, incubated at room temperature for 1hr.

Running condition: 10% gel, 130V for 1.5 hrs.


Transfer condition: Wet, 250mA, 2 hrs (NC membrane).

Blocking condition: 5% non-fat milk in TBS at 4°C overnight.

Washing condition: 5ml TBST, 4 x 5 min.

Exposure system: GE ECL 1 min.

Target bands were observed at 70kD in Human, mouse and rat cells and at 39kD in mouse cells.



Mouse kidney cells stained for ST2 using ab25877 at 20 µg/ml in ICC/IF.

Immunocytochemistry/ Immunofluorescence - Anti-ST2 antibody - BSA free (ab25877)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors