

SureLight® APC Anti-T7 tag® antibody ab72565

1 References [画像数 1](#)

製品の概要

製品名	SureLight® APC Anti-T7 tag® antibody
製品の詳細	SureLight® APC Rabbit polyclonal to T7 tag®
由来種	Rabbit
標識	SureLight® APC. Ex: 652nm, Em: 657nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Species independent
免疫原	Synthetic peptide T7 tag conjugated to a carrier molecule
特記事項	Excitation maximum:650 nm Emission maximum:660 nm The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	pH: 7.40 Preservative: 0.013% Sodium azide Constituents: 0.1% BSA, 0.87% Sodium chloride, 1.64% Sodium phosphate, 1.71% Sucrose
精製度	Size exclusion
特記事項(精製)	Size exclusion chromatography as well as some proprietary methods for purifying the final materials were used.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab72565の使用に適用されず
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent dilution.

ターゲット情報

関連性 The T7 tag is an 11 amino acid peptide encoded in the leader sequence of T7 bacteriophage gene10. This gene encodes a T7 major capsid protein whose function is not clear. The T7 tag serves as a tag in many expression vectors including the pET system that is based on the very efficient T7 RNA polymerase expression system. Monoclonal antibodies specific for T7 tag are an important tool for studying expression of recombinant T7-tagged proteins.

画像



Flow cytometry analysis of a recombinant (N-terminus T7-tag) transmembrane cell surface receptor in CHO cells using ab72565 at a 1/20 dilution (yellow). Black - Mouse IgG SureLight® Allophycocyanin isotype control.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors