

### Anti-SV40 T-antigen antibody [LP7N] ab43234

#### 製品の概要

製品名	Anti-SV40 T-antigen antibody [LP7N]
製品の詳細	Mouse monoclonal [LP7N] to SV40 T-antigen
由来種	Mouse
アプリケーション	<b>適用あり:</b> IHC-Fr
種交差性	<b>交差種:</b> Simian Virus 40
免疫原	Tissue, cells or virus corresponding to SV40 T-antigen. Detergent insoluble fraction of SVK14 cell line (SV40 transformed neonatal keratinocytes).
エピトープ	Epitope has been mapped to a region of the SV40 T antigen which is not covered by small t-shared domain or ND2 protein, since both are negative. Epitope is denaturation-sensitive and does not immunoblot.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
精製度	Protein A/G purified
ポリ/モノ	モノクローナル
クローン名	LP7N
アイソタイプ	IgG1

#### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab43234の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
IHC-Fr		Use at an assay dependent concentration.

## ターゲット情報

**細胞内局在**                      Nuclear

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors