


Anti-Synaptogyrin 1 antibody ab113886

1 References [画像数 2](#)

製品の概要

製品名	Anti-Synaptogyrin 1 antibody
製品の詳細	Rabbit polyclonal to Synaptogyrin 1
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Rat, Human 交差が予測される動物種: Mouse 
免疫原	Synthetic peptide from the internal region of Human Synaptogyrin 1, conjugated to KLH.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.4 Preservative: 0.02% Sodium azide Constituent: 99% PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、 次のテスト済みアプリケーションにおけるab113886の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 25 kDa.
IHC-P		Use a concentration of 5 µg/ml.

ターゲット情報

機能	Involved in the regulation of short-term and long-term synaptic plasticity.
配列類似性	Belongs to the synaptogyrin family. Contains 1 MARVEL domain.
細胞内局在	Membrane. Melanosome. Cell junction > synapse. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

画像



All lanes : Anti-Synaptogyrin 1 antibody (ab113886) at 1 µg/ml

Lane 1 : Rat brain tissue lysate

Lane 2 : Rat brain tissue lysate with with blocking peptide

Predicted band size: 25 kDa



ab113886 at 5ug/ml staining Synaptogyrin 1 in Human brain (cerebellum) tissue by immunohistochemistry (FFPE).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors