

Anti-Syndecan-1 antibody [B-A38], prediluted ab1410

2 References

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-Syndecan-1 antibody [B-A38], prediluted |
| 製品の詳細 | Mouse monoclonal [B-A38] to Syndecan-1, prediluted |
| 由来種 | Mouse |
| アプリケーション | 適用あり: IHC-Fr, IHC-P |
| 種交差性 | 交差種: Human |
| 免疫原 | Tissue, cells or virus corresponding to Syndecan-1. U266 cell line derived from Hodgkins disease |
| ポジティブ・コントロール | Tonsil |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | pH: 7.30 Constituents: PBS, Tissue culture supernatant |
| ポリモノ | モノクローナル |
| クローン名 | B-A38 |
| ミエローマ | unknown |
| アイソタイプ | IgG1 |
| 軽鎖の種類 | unknown |

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab1410の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------|
| IHC-Fr | | |
| IHC-P | | |

追加情報 IHC:
This antibody is pre-diluted for use with an ABC system.
For formalin fixed paraffin embedded sections, perform heat mediated antigen retrieval before commencing with IHC.

ターゲット情報

機能 Cell surface proteoglycan that bears both heparan sulfate and chondroitin sulfate and that links the cytoskeleton to the interstitial matrix.

配列類似性 Belongs to the syndecan proteoglycan family.

細胞内局在 Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors