


Anti-TACC2 antibody ab96760

面格粒

医薬用外劇物

製品の概要

製品名	Anti-TACC2 antibody
製品の詳細	Rabbit polyclonal to TACC2
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat 
免疫原	Synthetic peptide containing a sequence corresponding to a region within amino acids 237 and 303 of TACC2 (NP_996743; UniProt ID: O95359 isoform 5).
ポジティブ・コントロール	WB: H1299 whole cell lysate; 293T, A431, HeLaS3, HepG2, MOLT4 and Raji cell lysates IHC-P: Cal27 xenograft
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 1.21% Tris, 0.75% Glycine, 10% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab96760の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。


アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 119 kDa.
IHC-P		1/100 - 1/500.

ターゲット情報

関連性 Transforming acidic coiled coil proteins are a conserved family of centrosome and microtubule interacting proteins that are implicated in cancer. This gene encodes a protein that concentrates at centrosomes throughout the cell cycle. This gene lies within a chromosomal region associated with tumorigenesis. Expression of this gene is thought to affect the progression of breast tumors. Expression of this gene is also induced by erythropoietin.

細胞内局在 Cytoplasm. Nucleus. Cytoplasm, cytoskeleton, centrosome

画像



Western blot - Anti-TACC2 antibody (ab96760)

Anti-TACC2 antibody (ab96760) at 1/1000 dilution + H1299 whole cell lysate at 30 µg

Predicted band size: 119 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TACC2 antibody (ab96760)

Immunohistochemical analysis of TACC2 in paraffin embedded Cal27 xenograft, using ab96760 antibody at a 1/100 dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors