

Anti-TBX1 antibody [EPR3289(2)] ab109313

リコンビナント **RabMAb**

9 References **画像数 4**

製品の概要

製品名	Anti-TBX1 antibody [EPR3289(2)]
製品の詳細	Rabbit monoclonal [EPR3289(2)] to TBX1
由来種	Rabbit
アプリケーション	適用あり: WB 適用なし: Flow Cyt, ICC/IF or IHC-P
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: Jurkat, NCCIT, PC-12, F9, C6, Neuro-2a, HEK-293, HeLa cell lysates
特記事項	This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information see here . Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents .

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	pH: 7.20 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 40% Glycerol (glycerin, glycerine), 9.85% Tris glycine, 50% Tissue culture supernatant
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR3289(2)
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab109313の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/10000. Detects a band of approximately 43 kDa (predicted molecular weight: 43 kDa).

追加情報 Is unsuitable for Flow Cyt, ICC/IF or IHC-P.

ターゲット情報

機能	Probable transcriptional regulator involved in developmental processes. Is required for normal development of the pharyngeal arch arteries.
関連疾患	<p>Haploinsufficiency of the TBX1 gene is responsible for most of the physical malformations present in DiGeorge syndrome (DGS) and velocardiofacial syndrome (VCFS) [MIM:188400, 192430]. DGS is characterized by the association of several malformations: hypoplastic thymus and parathyroid glands, congenital conotruncal cardiopathy, and a subtle but characteristic facial dysmorphism. VCFS is marked by the association of congenital conotruncal heart defects, cleft palate or velar insufficiency, facial dysmorphology and learning difficulties. It is now accepted that these two syndromes represent two forms of clinical expression of the same entity manifesting at different stages of life.</p> <p>Defects in TBX1 are a cause of DiGeorge syndrome (DGS) [MIM:188400].</p> <p>Defects in TBX1 are a cause of velocardiofacial syndrome (VCFS) [MIM:192430].</p> <p>Defects in TBX1 are a cause of conotruncal heart malformations (CTHM) [MIM:217095]. CTHM consist of cardiac outflow tract defects, such as tetralogy of Fallot, pulmonary atresia, double-outlet right ventricle, truncus arteriosus communis, and aortic arch anomalies.</p>
配列類似性	Contains 1 T-box DNA-binding domain.
細胞内局在	Nucleus.

画像



All lanes : Anti-TBX1 antibody [EPR3289(2)] (ab109313) at 1/1000 dilution

Lane 1 : HEK-293 (Human embryonic kidney epithelial cell) control cell lysate

Lane 2 : HEK-293 (Human embryonic kidney epithelial cell) nuclear cell lysate

Lane 3 : HeLa (Human cervix adenocarcinoma epithelial cell) whole cell lysate

Lane 4 : Human liver lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG (HRP) with minimal cross-reactivity with human IgG at 1/20000 dilution


Predicted band size: 43 kDa

Observed band size: 53 kDa

Exposure time: 10 seconds

Blocking buffer: 5% NFDM/TBST.

There is a strong non-specific band in human liver. Human liver is negative control for TBX1 which is consistent with PMID: 28920943.



All lanes : Anti-TBX1 antibody [EPR3289(2)] (ab109313) at 1/2000 dilution

Lane 1 : Jurkat (Human T cell leukemia T lymphocyte) whole cell lysate

Lane 2 : NCCIT (Human pluripotent embryonic carcinoma epithelial cell) whole cell lysate

Lane 3 : PC-12 (Rat adrenal gland pheochromocytoma) whole cell lysate

Lane 4 : F9 (Mouse embryonal carcinoma epithelial cell) whole cell lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000


dilution

Predicted band size: 43 kDa

Observed band size: 53 kDa

Exposure time: 30 seconds

Blocking buffer: 5% NFDM/TBST.



Lane 1 : Anti-TBX1 antibody [EPR3289(2)] (ab109313) at 1/2000 dilution

Lane 2 : Anti-TBX1 antibody [EPR3289(2)] (ab109313) at 1/5000 dilution

Lane 1 : C6 (Rat glial tumor glial cell) whole cell lysate

Lane 2 : Neuro-2a (Mouse neuroblastoma neuroblast) whole cell lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution

Predicted band size: 43 kDa

Observed band size: 53 kDa

Blocking buffer: 5% NFDM/TBST.

Exposure time: 30 seconds (Lane 1) and 180 seconds (Lane 2).

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-TBX1 antibody [EPR3289(2)] (ab109313)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors