

Anti-Testosterone 3 CMO antibody [3T16] ab20250

製品の概要

製品名	Anti-Testosterone 3 CMO antibody [3T16]
製品の詳細	Mouse monoclonal [3T16] to Testosterone 3 CMO
由来種	Mouse
アプリケーション	適用あり: RIA
種交差性	交差種: Species independent
免疫原	Chemical/ Small Molecule corresponding to Testosterone 3 CMO conjugated to bovine serum albumin. Testosterone 3 CMO conjugated to BSA.
特記事項	<p>This antibody can be used for detecting and measuring testosterone levels in serum and tissue.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
精製度	Protein A/G purified
一次抗体 備考	This antibody can be used for detecting and measuring testosterone levels in serum and tissue.
ポリ/モノ	モノクローナル
クローン名	3T16
ミエローマ	P3-NS1/1-Ag4-1
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab20250の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
RIA		Use at an assay dependent dilution.

ターゲット情報

関連性 Testosterone-3-carboxy-methyl-oxyme (T-3-CMO). Testosterone belongs to a class of hormones called androgens. These are male hormones responsible for the development of the male reproductive system and secondary male sexual characteristics such as voice depth and facial hair. Testosterone is normally produced by the testes in large quantities in men. It also occurs normally in smaller quantities in women. Several man-made derivatives of testosterone are used to treat advanced disseminated breast cancer in women, especially when cancer has spread to the bones.

細胞内局在 Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors