

# Anti-Tetanus Toxoid antibody ab24292

### 製品の概要

製品名	Anti-Tetanus Toxoid antibody
製品の詳細	Goat polyclonal to Tetanus Toxoid
由来種	Goat
アプリケーション	<b>適用あり:</b> Immunoelectrophoresis
種交差性	<b>交差種:</b> Clostridium tetani
免疫原	Tetanus toxoid (purified)
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.1% Sodium azide Constituent: 0.0268% PBS
精製度	Ion Exchange Chromatography
特記事項 (精製)	>95% pure. Sodium sulfate precipitation and ion-exchange chromatography
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab24292の使用に適用されます  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Immunoelectrophoresis		

#### 追加情報

le: Use at an assay dependent dilution.  
May be used in place of neat antiserum in almost any appropriate antibody-based technique.

Not tested in other applications.  
Optimal dilutions/concentrations should be determined by the end user.

#### ターゲット情報

#### 関連性

Tetanus toxoid is derived from the toxin released by Clostridium Tetanis that causes the disease tetanus. It is used as a vaccine to prevent tetanus or to help boost the immune response to other vaccines. A toxoid is a substance that has been treated to destroy its toxic properties but that retains the capacity to stimulate production of antitoxins, used in immunization. Tetanus toxoid is produced by formaldehyde inactivation of pure neuro toxin. Each lot is verified to be non toxic. In addition to being a potent antigen, tetanus toxoid is a useful carrier protein for making polysaccharides and haptens immunogenic.

#### 細胞内局在

Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors