

### Anti-TGE Virus peplomer antibody [1Q17] ab2599

#### 製品の概要

---

製品名	Anti-TGE Virus peplomer antibody [1Q17]
製品の詳細	Mouse monoclonal [1Q17] to TGE Virus peplomer
由来種	Mouse
特異性	This antibody is specific for TGE (Transmissible gastroenteritis) virus peplomer antigen.
アプリケーション	<b>適用あり:</b> WB, ELISA
種交差性	<b>交差種:</b> Transmissible gastroenteritis
免疫原	Tissue, cells or virus corresponding to TGE Virus peplomer. Purified TGE virus
エピトープ	The exact epitope recognised by this antibody has not been mapped.
特記事項	<p>The peplomer protein mediates the binding of virions to the host cell receptor and is also involved with membrane fusion.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.20 Preservative: 0.01% Sodium azide
一次抗体 備考	The peplomer protein mediates the binding of virions to the host cell receptor and is also involved with membrane fusion.
ポリ/モノ	モノクローナル
クローン名	1Q17
ミエローマ	unknown

アイソタイプ	IgG1
軽鎖の種類	unknown

## アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab2599の使用に適用されます**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/1000.
ELISA		1/10000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors