

Anti-TGF alpha antibody [P/T1] ab9578

★★★★☆ [1 Abreviews](#) [6 References](#) [画像数 2](#)

製品の概要

製品名	Anti-TGF alpha antibody [P/T1]
製品の詳細	Mouse monoclonal [P/T1] to TGF alpha
由来種	Mouse
特異性	This antibody was selected and evaluated immunoserologically for recombinant human TGF-alpha and synthetic (uncoupled) peptides using simple binding and competitive ELISA. Immunohistochemical characterization of the anti-TGF-alpha monoclonal antibody was made on different human carcinomas. Recommended for immunohistochemical detection of TGF-alpha production in human tumor cells, immunoserological monitoring of progression of malignant tumours and for screening of therapeutic efficiency of different human tumors.
アプリケーション	適用あり: IHC-P, Flow Cyt
種交差性	交差種: Human
免疫原	Synthetic peptide fragment (32-40 amino acids sequence of 50 amino acids) of transforming growth factor alpha (TGF-alpha).
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.097% Sodium azide Constituents: Tris HCl, 0.2% BSA
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	P/T1

アイソタイプ IgG1
軽鎖の種類 kappa

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab9578の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P	★★★★★ (1)	1/100. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.
Flow Cyt		Use 1µg for 10 ⁶ cells. ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

ターゲット情報


機能 TGF alpha is a mitogenic polypeptide that is able to bind to the EGF receptor and to act synergistically with TGF beta to promote anchorage-independent cell proliferation in soft agar.

組織特異性 Isoform 1, isoform 3 and isoform 4 are expressed in keratinocytes and tumor-derived cell lines.

配列類似性 Contains 1 EGF-like domain.


細胞内局在 Cell membrane and Secreted > extracellular space.

画像



Immunohistochemical characterization of the anti-TGF-a monoclonal antibody ab9578 on human melanoma.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TGF alpha antibody [P/T1] (ab9578)



Overlay histogram showing A431 cells stained with ab9578 (red line). The cells were fixed with 80% methanol (5 min) and incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions. The cells were then incubated with the antibody (ab9578, 1µg/1x10⁶ cells) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-mouse IgG (H+L) ([ab96879](#)) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG1 [ICIGG1] ([ab91353](#), 2µg/1x10⁶ cells) used under the same conditions. Acquisition of >5,000 events was performed.

Please note that Abcam do not have any data for use of this antibody on non-fixed cells. We welcome any customer feedback.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors