

### Anti-TGM1 antibody ab27000

★☆☆☆☆ [2 Abreviews](#) [3 References](#) [画像数 1](#)

#### 製品の概要

製品名	Anti-TGM1 antibody
製品の詳細	Rabbit polyclonal to TGM1
由来種	Rabbit
特異性	This antibody is specific to keratinocyte Transglutaminase 1.
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide corresponding to Human TGM1. A 15 residue synthetic peptide derived from a domain specific for TGM1. Database link: <a href="#">P22735</a>
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Constituent: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

The Abpromise guarantee

**Abpromise保証は、次のテスト済みアプリケーションにおけるab27000の使用に適用されます**


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB	★☆☆☆☆ (1)	1/200 - 1/1000. Predicted molecular weight: 90 kDa.

## ターゲット情報

機能	Catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins. Responsible for cross-linking epidermal proteins during formation of the stratum corneum.
関連疾患	<p>Defects in TGM1 are the cause of ichthyosis lamellar type 1 (LI1) [MIM:242300]. LI is a non-bullous ichthyosis, a skin disorder characterized by abnormal cornification of the epidermis. It is one the most severe forms of ichthyoses apparent at birth and persisting throughout life. LI patients are born encased in a tight, shiny, translucent covering called collodion membrane. Over the first weeks of life, the collodion membrane is gradually replaced by generalized large, dark brown, plate-like scales with minimal to no erythroderma. Tautness of facial skin commonly results in ectropion, eclabium and scarring alopecia of the scalp. Common complications are severe heat intolerance and recurrent ear infections.</p> <p>Defects in TGM1 are a cause of non-bullous congenital ichthyosiform erythroderma (NCIE) [MIM:242100]. NCIE is a non-bullous ichthyosis, a skin disorder characterized by abnormal cornification of the epidermis. Most affected individuals are born with a tight, shiny, translucent covering called collodion membrane. The collodion membrane subsequently evolves into generalized scaling and intense redness of the skin. Clinical features are milder than in lamellar ichthyoses and demonstrate a greater variability in the intensity of erythema, size and type of scales. In contrast to lamellar ichthyoses, scales are usually white, fine and powdery, and palms and soles are severely affected. Patients suffer from palmoplantar keratoderma, often with painful fissures, digital contractures, and loss of pulp volume.</p> <p>Defects in TGM1 are the cause of ichthyosis congenital autosomal recessive TGM1-related (ARCI-TGM1) [MIM:242300]. A disorder of keratinization with abnormal differentiation and desquamation of the epidermis resulting in two major clinical entities. Lamellar ichthyosis is a condition often associated with an embedment in a collodion-like membrane at birth; skin scales later develop, covering the entire body surface. Non-bullous congenital ichthyosiform erythroderma characterized by fine whitish scaling on an erythrodermal background; larger brownish scales are present on the buttocks, neck and legs.</p>
配列類似性	Belongs to the transglutaminase superfamily. Transglutaminase family.
翻訳後修飾	The membrane anchorage region possesses a cluster of five cysteines within which fatty acid(s) may become thioester-linked. It is subject to phorbol ester-stimulated phosphorylation and is hypersensitive to proteolysis, which releases the enzyme in a soluble form.
細胞内局在	Membrane.

## 画像



Anti-TGM1 antibody (ab27000) at 1/500 dilution + Human recombinant Transglutaminase at 0.5 µg

**Secondary**

HRP-conjugated anti Rabbit IgG at 1/5000 dilution

**Predicted band size:** 90 kDa

**Observed band size:** 90 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors