


# Anti-Thrombospondin 2 antibody ab84469

**11 References**   [画像数 2](#)

### 製品の概要

製品名	Anti-Thrombospondin 2 antibody
製品の詳細	Rabbit polyclonal to Thrombospondin 2
由来種	Rabbit
特異性	Reactivity in Rat is no longer guaranteed due to difficulty reproducing in-house experimental evidence, but will remain predicted, due to the immunogen showing 100% xr with Rat. From Jan 2024, QC testing of replenishment batches of this polyclonal changed. All tested and expected application and reactive species combinations are still covered by our Abcam product promise. However, we no longer test all applications. For more information on a specific batch, please contact our Scientific Support who will be happy to help.
アプリケーション	<b>適用あり:</b> WB, ICC/IF
種交差性	<b>交差種:</b> Mouse, Human <b>交差が予測される動物種:</b> Rat, Horse, Chicken, Cow, Dog 
免疫原	Synthetic peptide corresponding to Human Thrombospondin 2 aa 650-750 conjugated to keyhole limpet haemocyanin. (Peptide available as <a href="#">ab96113</a> )
ポジティブ・コントロール	This antibody gave a positive signal in both Human and Mouse spinal cord tissue lysate. This antibody gave a positive signal in ICC/IF in U2O2 cells.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40

Preservative: 0.02% Sodium azide  
Constituent: PBS

Batches of this product that have a concentration < 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our scientific support team who will be happy to help.

精製度 Immunogen affinity purified  
ポリ/モノ ポリクローナル  
アイソタイプ IgG

## アプリケーション

**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab84469の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 58 kDa (predicted molecular weight: 129 kDa).
ICC/IF		Use a concentration of 10 µg/ml.

## ターゲット情報


**機能** Adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. Can bind to fibrinogen, fibronectin, laminin and type V collagen.

**組織特異性** High expression in intervertebral disk tissue.

**関連疾患** Genetic variations in THBS2 may be a cause of susceptibility to intervertebral disc disease (IDD) [MIM:603932]; also known as lumbar disk herniation (LDH). IDD is one of the most common musculo-skeletal disorders and the predominant cause of low-back pain and unilateral leg pain.


**配列類似性** Belongs to the thrombospondin family.  
Contains 3 EGF-like domains.  
Contains 1 TSP C-terminal (TSPC) domain.  
Contains 1 TSP N-terminal (TSPN) domain.  
Contains 3 TSP type-1 domains.  
Contains 8 TSP type-3 repeats.  
Contains 1 VWFC domain.

## 画像



Immunocytochemistry/ Immunofluorescence - Anti-Thrombospondin 2 antibody (ab84469)

Ab84469 staining Thrombospondin 2 in U2O2 cells. The cells were fixed with 100% Methanol (5 min), permeabilized with 0.1% Triton X-100 for 5 minutes and then blocked with 1% BSA/10% normal goat serum/0.3M glycine in 0.1% PBS-Tween for 1h. The cells were then incubated overnight at +4°C with ab84469 at 10ug/ml then detected with an Alexa Fluor® 488 goat anti-rabbit secondary antibody (**ab150081**) at a 1/1000 dilution (shown in green). Nuclear DNA was labelled with DAPI (shown in blue), and tubulin was labelled using **ab7291** Mouse monoclonal to alpha Tubulin together with **ab150120** Goat Anti-Mouse IgG H&L (Alexa Fluor® 594), at a 1/1000 dilution (shown in red).



Western blot - Anti-Thrombospondin 2 antibody (ab84469)

**All lanes :** Anti-Thrombospondin 2 antibody (ab84469) at 1 mg/ml

**Lane 1 :** Mouse spinal cord embryo E14

**Lane 2 :** Human skin

**Lane 3 :** Mouse skin

**Lane 4 :** A431 (Human epithelial carcinoma cell line)

Lysates/proteins at 10 µg per lane.

### Secondary

**All lanes :** Goat polyclonal to Rabbit IgG - H&L - Pre-Adsorbed (HRP) at 1/50000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

**Predicted band size:** 129 kDa

**Additional bands at:** 210 kDa (possible glycosylated form), 46 kDa (possible cleavage fragment), 58 kDa (possible non-specific binding)

**Exposure time:** 20 minutes

This blot was produced using a 3-8% Tris Acetate gel under the TA buffer system. The gel was run at 150V for 60 minutes before being transferred onto a Nitrocellulose membrane at 30V for 70 minutes. The membrane was then blocked for an hour using 2% Bovine Serum Albumin before being incubated with abX overnight at 4°C. Antibody binding was detected using an anti-rabbit antibody conjugated to HRP, and visualised using ECL development solution **ab133406**.

ab84469 contains a potential glycosylation site (SwissProt) which may explain its migration at a higher molecular weight than predicted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors