

Anti-Thymidylate Synthase antibody [EPR4545] ab108995

リコンビナント **RabMAb**

★★★★★ **2 Abreviews** **18 References** 画像数 8

製品の概要

製品名	Anti-Thymidylate Synthase antibody [EPR4545]
製品の詳細	Rabbit monoclonal [EPR4545] to Thymidylate Synthase
由来種	Rabbit
特異性	The mouse and rat recommendation is based on the WB results. We do not guarantee IHC-P for mouse and rat.
アプリケーション	適用あり: Flow Cyt (Intra), WB, IP, IHC-P, ICC/IF
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	HeLa, Jurkat, HL60 and MCF7 cell lysates; Human ovary carcinoma and colon carcinoma tissue.
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.20</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 40% Glycerol, PBS, 0.05% BSA</p>
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR4545

アプリケーション


The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab108995の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt (Intra)		1/10 - 1/100. ab172730 - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.
WB		1/1000 - 1/10000. Predicted molecular weight: 36 kDa.
IP		1/10 - 1/100.
IHC-P	★★★★★ (1)	1/50 - 1/100. Perform heat mediated antigen retrieval before commencing with IHC staining protocol. Heat up to 98°C, below boiling, and then let cool for 10-20 minutes. The mouse and rat recommendation is based on the WB results. We do not guarantee IHC-P for mouse and rat.
ICC/IF		1/100 - 1/250.

ターゲット情報

機能	Contributes to the de novo mitochondrial thymidylate biosynthesis pathway.
パスウェイ	Pyrimidine metabolism; dTTP biosynthesis.
配列類似性	Belongs to the thymidylate synthase family.
細胞内局在	Nucleus. Cytoplasm. Mitochondrion. Mitochondrion matrix. Mitochondrion inner membrane.

画像



All lanes : Anti-Thymidylate Synthase antibody [EPR4545] (ab108995) at 1/10000 dilution (purified)

Lane 1 : HeLa (Human cervix adenocarcinoma epithelial cell) whole cell lysates

Lane 2 : Jurkat (Human T cell leukemia T lymphocyte) whole cell lysates

Lane 3 : C6 (Rat glial tumor glial cell) whole cell lysates

Lane 4 : RAW 264.7 (Mouse Abelson murine leukemia virus-induced tumor macrophage) whole cell lysates

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000 dilution

Predicted band size: 36 kDa

Blocking and diluting buffer: 5% NFDM/TBST.

ab108995 (purified) at 1:30 dilution (2µg) immunoprecipitating Thymidylate Synthase in Jurkat whole cell lysate.


Lane 1 (input): Jurkat (Human T cell leukemia T lymphocyte) whole cell lysate, 10µg


Lane 2 (+): ab108995 & Jurkat whole cell lysate

Lane 3 (-): Rabbit monoclonal IgG (**ab172730**) instead of ab108995 in Jurkat whole cell lysate

For western blotting, VeriBlot for IP Detection Reagent (HRP) (**ab131366**) was used for detection at 1:1000 dilution.


Blocking and diluting buffer: 5% NFDM/TBST.






Flow Cytometry (Intracellular) - Anti-Thymidylate Synthase antibody [EPR4545] (ab108995)

Intracellular Flow Cytometry analysis of HeLa (Human cervix adenocarcinoma epithelial cell) cells labeling Thymidylate Synthase with purified ab108995 at 1/50 dilution (red). Cells were fixed with 4% Paraformaldehyde. A Goat anti rabbit IgG (Alexa Fluor® 488) secondary antibody was used at 1/2000 dilution. Isotype control - Rabbit monoclonal IgG (Black). Unlabeled control - Cell without incubation with primary antibody and secondary antibody (Blue).




Immunocytochemistry/ Immunofluorescence - Anti-Thymidylate Synthase antibody [EPR4545] (ab108995)

Immunocytochemistry/ Immunofluorescence analysis of Jurkat (Human T cell leukemia T lymphocyte) cells labeling Thymidylate Synthase with Purified ab108995 at 1:100 dilution. Cells were fixed in 4% Paraformaldehyde and permeabilized with 0.1% tritonX-100 **ab150077** Goat anti rabbit IgG(Alexa Fluor® 488) was used as the secondary antibody at 1:1000 dilution. DAPI nuclear counterstain. PBS instead of the primary antibody was used as the secondary antibody only control.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Thymidylate Synthase antibody [EPR4545] (ab108995)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human colon carcinoma tissue sections labeling Thymidylate Synthase with purified ab108995 at 1:100 dilution. Heat mediated antigen retrieval was performed using **ab93684** (Tris/EDTA buffer, pH 9.0). Tissue was counterstained with hematoxylin. ImmunoHistoProbe one step HRP Polymer (ready to use) secondary antibody was used at 1:0 dilution. PBS instead of the primary antibody was used as the negative control.



Western blot - Anti-Thymidylate Synthase antibody [EPR4545] (ab108995)

All lanes : Anti-Thymidylate Synthase antibody [EPR4545] (ab108995) at 1/1000 dilution (unpurified)

Lane 1 : HeLa cell lysate

Lane 2 : Jurkat cell lysate


Lane 3 : HL60 cell lysate

Lane 4 : MCF7 cell lysate

Lysates/proteins at 10 µg per lane.

Performed under reducing conditions.

Predicted band size: 36 kDa




Unpurified ab108995 at 1/50 dilution, staining Thymidylate Synthase in Human ovary carcinoma Immunohistochemistry, Paraffin-embedded tissue.

Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Thymidylate Synthase antibody [EPR4545] (ab108995)

Why choose a recombinant antibody?



Anti-Thymidylate Synthase antibody [EPR4545] (ab108995)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors