

### Anti-TMEM184A antibody ab106702

画像数 1

#### 製品の概要

|              |   |
|--------------|---|
| 製品名          | Anti-TMEM184A antibody  |
| 製品の詳細        | Rabbit polyclonal to TMEM184A   |
| 由来種          | Rabbit  |
| 特異性          | ab106702 is predicted to not cross-react with other members of the TMEM184 family.  |
| アプリケーション     | <b>適用あり:</b> WB   |
| 種交差性         | <b>交差種:</b> Rat   |
| 免疫原          | A 16 amino acid synthetis peptide near the carboxy terminus of Human TMEM184A (NP_001091089)  |
| ポジティブ・コントロール | Rat brain tissue lysate   |
| 特記事項         | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

#### 製品の特性

|        |  |
|--------|--|
| 製品の状態  | Liquid   |
| 保存方法   | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle. Stable for 12 months at -20°C. |
| バッファー  | pH: 7.2<br>Preservative: 0.02% Sodium azide<br>Constituent: PBS  |
| 精製度    | Immunogen affinity purified  |
| ポリ/モノ  | ポリクローナル  |
| アイソタイプ | IgG  |

#### アプリケーション


**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab106702の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

| アプリケーション | Abreviews | 特記事項  |
|----------|-----------|---|
| WB       |           | Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 46 kDa. |

#### ターゲット情報

**配列類似性**                      Belongs to the TMEM184 family.  
**細胞内局在**                      Membrane.

#### 画像



Western Blot of rat brain tissue lysate labeling TMEM184A with Anti-TMEM184A antibody (ab106702) at 1µg/ml in (1) the absence and (2) the presence of blocking peptide.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors