

Anti-TNFAIP3 antibody [EPR2663] ab92324

KO 評価済 リコンビナント RabMAb

★★★★★ 1 Abreviews 31 References 画像数 10

製品の概要

製品名	Anti-TNFAIP3 antibody [EPR2663]
製品の詳細	Rabbit monoclonal [EPR2663] to TNFAIP3
由来種	Rabbit
特異性	Mouse species is recommended based on WB results, we do not guarantee IHC-P for mouse.
アプリケーション	適用あり: WB, IHC-P 適用なし: Flow Cyt (Intra) or ICC/IF
種交差性	交差種: Mouse, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: WEHI-3 treated with TNF (ab9642), Jurkat treated with TNF (ab9642) + TPA, Jurkat treated with 5ng/ml PMA for 48 hours and then treated with 2 μ g/ml PHA for 48 hours, HeLa, A549 and Daudi cell lysates. IHC-P: Human kidney tissue.
特記事項	This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information see here . Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents .

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Stable for 12 months at -20°C.
バッファー	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 59% PBS, 40% Glycerol, 0.05% BSA
精製度	Protein A purified

ポリクローナル	モノクローナル
クローン名	EPR2663
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab92324の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/5000. Predicted molecular weight: 90 kDa.
IHC-P	★★★★★ (1)	1/50 - 1/100. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol. Mouse species is recommended based on WB results, we do not guarantee IHC-P for mouse.

追加情報 Is unsuitable for Flow Cyt (Intra) or ICC/IF.

ターゲット情報

機能	Ubiquitin-editing enzyme that contains both ubiquitin ligase and deubiquitinase activities. Essential component of a ubiquitin-editing protein complex, comprising also RNF11, ITCH and TAX1BP1, that ensures the transient nature of inflammatory signaling pathways. Upon TNF stimulation, deubiquitinates 'Lys-63'-polyubiquitin chains on RIPK1 and catalyzes the formation of 'Lys-48'-polyubiquitin chains. This leads to RIPK1 proteasomal degradation and consequently termination of the TNF- or LPS-mediated activation of NF-kappa-B. In vitro able to deubiquitinate both 'Lys-48'- and 'Lys-63' polyubiquitin chains. Inhibitor of programmed cell death. Has a role in the function of the lymphoid system.
配列類似性	Belongs to the peptidase C64 family. Contains 7 A20-type zinc fingers. Contains 1 OTU domain.
ドメイン	The A20-type zinc fingers mediate the ubiquitin ligase activity. The OTU domain mediates the deubiquitinase activity.
細胞内局在	Cytoplasm. Nucleus.

画像



Western blot - Anti-TNFAIP3 antibody [EPR2663]
(ab92324)

All lanes : Anti-TNFAIP3 antibody [EPR2663] (ab92324) at 1/1000 dilution

Lane 1 : Wild-type A549 cell lysate

Lane 2 : TNFAIP3 knockout A549 cell lysate

Lane 3 : Wild-type HeLa cell lysate

Lane 4 : TNFAIP3 knockout HeLa cell lysate

Lysates/proteins at 20 µg per lane.


Performed under reducing conditions.

Predicted band size: 90 kDa

Observed band size: 90 kDa


Lanes 1- 4: Merged signal (red and green). Green - ab92324 observed at 90 kDa. Red - Anti-GAPDH antibody [6C5] - Loading Control ([ab8245](#)) observed at 37 kDa.

ab92324 was shown to react with TNFAIP3 in wild-type A549 cells in western blot. Loss of signal was observed when knockout cell line [ab266946](#) (knockout cell lysate [ab257114](#)) was used. Wild-type A549 and TNFAIP3 knockout A549 cell lysates were subjected to SDS-PAGE. Membrane was blocked for 1 hour at room temperature in 0.1% TBST with 3% non-fat dried milk. ab92324 and Anti-GAPDH antibody [6C5] - Loading Control ([ab8245](#)) overnight at 4°C at a 1 in 1000 dilution and a 1 in 20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye®800CW) preadsorbed ([ab216773](#)) and Goat anti-Mouse IgG H&L (IRDye®680RD) preadsorbed ([ab216776](#)) secondary antibodies at 1 in 20000 dilution for 1 hour at room temperature before imaging.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TNFAIP3 antibody [EPR2663] (ab92324)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human kidney tissue labelling TNFAIP3 with unpurified ab92324 at 1/50. Heat mediated antigen retrieval was performed using Tris/EDTA buffer pH 9. A prediluted HRP-polymer conjugated anti-rabbit IgG was used as the secondary antibody. Counterstained with Hematoxylin.



Western blot - Anti-TNFAIP3 antibody [EPR2663] (ab92324)

All lanes : Anti-TNFAIP3 antibody [EPR2663] (ab92324) at 1/1000 dilution

Lane 1 : Wild-type HeLa cell lysate

Lane 2 : TNFAIP3 knockout HeLa cell lysate

Lane 3 : Jurkat cell treated with 5ng/ml PMA for 48 hours and then treated with 2µg/ml PHA for 48 hours, whole cell lysate

Lane 4 : Untreated Jurkat cell lysate

Lysates/proteins at 20 µg per lane.

Performed under reducing conditions.


Predicted band size: 90 kDa

Observed band size: 80 kDa

Lanes 1-4: Merged signal (red and green). Green - ab92324 observed at 80 kDa. Red - loading control, **ab8245** observed at 37 kDa.

ab92324 Anti-TNFAIP3 antibody [EPR2663] was shown to specifically react with TNFAIP3 in wild-type HeLa cells. Loss of signal was observed when knockout cell line **ab265983** (knockout

cell lysate [ab257112](#)) was used. Wild-type and TNFAIP3 knockout samples were subjected to SDS-PAGE. ab92324 and Anti-GAPDH antibody [6C5] - Loading Control ([ab8245](#)) were incubated overnight at 4°C at 1 in 1000 dilution and 1 in 20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed ([ab216773](#)) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed ([ab216776](#)) secondary antibodies at 1 in 10000 dilution for 1 hour at room temperature before imaging.



Western blot - Anti-TNFAIP3 antibody [EPR2663] (ab92324)

All lanes : Anti-TNFAIP3 antibody [EPR2663] (ab92324) at 1/1000 dilution

Lane 1 : Wild-type A549 (Human lung carcinoma cell line) whole cell lysate

Lane 2 : TNFAIP3 knockout A549 (Human lung carcinoma cell line) whole cell lysate

Lane 3 : Jurkat (Human T cell leukemia cell line from peripheral blood) cell treated with 5ng/ml PMA for 48 hours and then treated with 2µg/ml PHA for 48 hours, whole cell lysate

Lane 4 : Untreated Jurkat (Human T cell leukemia cell line from peripheral blood) whole cell lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed ([ab216773](#)) at 1/10000 dilution


Predicted band size: 90 kDa

Observed band size: 80 kDa

Lanes 1-4: Merged signal (red and green). Green - ab92324 observed at 80 kDa. Red - loading control [ab8245](#) observed at 36 kDa.

ab92324 Anti-TNFAIP3 antibody [EPR2663] was shown to specifically react with TNFAIP3 in wild-type A549 cells. Loss of signal was observed when knockout cell line [ab266945](#) (knockout cell lysate [ab257113](#)) was used. Wild-type and TNFAIP3 knockout samples were subjected to SDS-PAGE. ab92324 and Anti-GAPDH antibody [6C5] - Loading Control ([ab8245](#)) were incubated

overnight at 4°C at 1 in 1000 dilution and 1 in 20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed ([ab216773](#)) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed ([ab216776](#)) secondary antibodies at 1 in 20000 dilution for 1 hour at room temperature before imaging.



Western blot - Anti-TNFAIP3 antibody [EPR2663] ([ab92324](#))

All lanes : Anti-TNFAIP3 antibody [EPR2663] ([ab92324](#)) at 1/5000 dilution

Lane 1 : WEHI-3 (Mouse leukemia lymphoblast) whole cell lysate

Lane 2 : WEHI-3 treated with 20 ng/ml TNF alpha ([ab9642](#)) for 6 h


Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution

Predicted band size: 90 kDa

Observed band size: 80 kDa



Western blot - Anti-TNFAIP3 antibody [EPR2663] ([ab92324](#))

Anti-TNFAIP3 antibody [EPR2663] ([ab92324](#)) at 1/2000 dilution (unpurified) + Jurkat cell lysate - treated with TNF ([ab9642](#)) and TPA at 10 µg


Secondary

Peroxidase-conjugated goat anti-rabbit IgG (H+L) at 1/1000 dilution

Predicted band size: 90 kDa

Observed band size: 80 kDa

Blocking and dilution buffer: 5% NFDm/TBST.



Western blot - Anti-TNFAIP3 antibody [EPR2663] (ab92324)

Anti-TNFAIP3 antibody [EPR2663] (ab92324) at 1/3000 dilution (purified) + Jurkat cell lysate - treated with TNF ([ab9642](#)) and TPA at 10 µg


Secondary

Peroxidase-conjugated goat anti-rabbit IgG (H+L) at 1/1000 dilution

Predicted band size: 90 kDa


Observed band size: 80 kDa

Blocking and dilution buffer: 5% NFDM/TBST.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TNFAIP3 antibody [EPR2663] (ab92324)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human kidney tissue labelling TNFAIP3 with purified ab92324 at 1/100. Heat mediated antigen retrieval was performed using Tris/EDTA buffer pH 9. A prediluted HRP-polymer conjugated anti-rabbit IgG was used as the secondary antibody. Counterstained with Hematoxylin.



Western blot - Anti-TNFAIP3 antibody [EPR2663] (ab92324)

All lanes : Anti-TNFAIP3 antibody [EPR2663] (ab92324) at 1/1000 dilution (unpurified)

Lane 1 : Jurkat cells treated with TNF ([ab9642](#)) and TPA

Lane 2 : Daudi cell lysate

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : HRP labelled goat anti-rabbit IgG at 1/2000 dilution

Predicted band size: 90 kDa

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-TNFAIP3 antibody [EPR2663] (ab92324)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors