

# Anti-Transferrin Receptor 2/TFR2 antibody ab80194

★★★★★ [32 Abreviews](#) [5 References](#) [画像数 3](#)

### 製品の概要

製品名	Anti-Transferrin Receptor 2/TFR2 antibody
製品の詳細	Rabbit polyclonal to Transferrin Receptor 2/TFR2
由来種	Rabbit
特異性	From Jan 2024, QC testing of replenishment batches of this polyclonal changed. All tested and expected application and reactive species combinations are still covered by our Abcam product promise. However, we no longer test all applications. For more information on a specific batch, please contact our Scientific Support who will be happy to help.
アプリケーション	<b>適用あり:</b> IHC-P, IHC-Fr, WB <b>適用なし:</b> ICC/IF
種交差性	<b>交差種:</b> Mouse, Rat, Human
免疫原	Synthetic peptide within Human Transferrin Receptor 2/TFR2 aa 150-250 conjugated to keyhole limpet haemocyanin. The exact sequence is proprietary. (Peptide available as <a href="#">ab87613</a> )
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40 Preservative: 0.02% Sodium azide Constituent: PBS
	Batches of this product that have a concentration < 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our

	scientific support team who will be happy to help.
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション

**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab80194の使用に適用されます  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 5 - 10 µg/ml. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.
IHC-Fr	★★★★★ (10)	Use at an assay dependent concentration.
WB	★★★★★ (20)	Use a concentration of 1 µg/ml. Detects a band of approximately 88 kDa (predicted molecular weight: 88 kDa). Abcam recommends using 5% BSA as the blocking agent.

**追加情報** Is unsuitable for ICC/IF.

## ターゲット情報

機能	Mediates cellular uptake of transferrin-bound iron in a non-iron dependent manner. May be involved in iron metabolism, hepatocyte function and erythrocyte differentiation.
組織特異性	Predominantly expressed in liver. While the alpha form is also expressed in spleen, lung, muscle, prostate and peripheral blood mononuclear cells, the beta form is expressed in all tissues tested, albeit weakly.
関連疾患	Defects in TFR2 are a cause of hemochromatosis type 3 (HFE3) [MIM:604250]. HFE3 is a disorder of iron hemostasis resulting in iron overload and has a phenotype indistinguishable from that of hereditary hemochromatosis (HH). HH is characterized by abnormal intestinal iron absorption and progressive increase of total body iron, which results in midlife in clinical complications including cirrhosis, cardiopathy, diabetes, endocrine dysfunctions, arthropathy, and susceptibility to liver cancer. Since the disease complications can be effectively prevented by regular phlebotomies, early diagnosis is most important to provide a normal life expectancy to the affected subjects.
配列類似性	Belongs to the peptidase M28 family. M28B subfamily.
細胞内局在	Cell membrane and Cytoplasm. Lacks the transmembrane domain. Probably intracellular.

## 画像



Western blot - Anti-Transferrin Receptor 2/TFR2 antibody (ab80194)

Anti-Transferrin Receptor 2/TFR2 antibody (ab80194) at 1 µg/ml + Human Liver Tissue Lysate at 20 µg

### Secondary

Goat polyclonal to Rabbit IgG - H&L - Pre-Adsorbed (HRP) at 1/3000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

**Predicted band size:** 88 kDa

**Observed band size:** 89 kDa

**Additional bands at:** 50 kDa, 98 kDa (possible glycosylated form). We are unsure as to the identity of these extra bands.

**Exposure time:** 8 minutes



Western blot - Anti-Transferrin Receptor 2/TFR2 antibody (ab80194)

This image is courtesy of an anonymous Abreview

**All lanes :** Anti-Transferrin Receptor 2/TFR2 antibody (ab80194) at 1/1000 dilution

**All lanes :** Rat liver whole tissue lysate

### Secondary

**All lanes :** HRP-conjugated goat anti-rabbit IgG polyclonal at 1/2000 dilution


Developed using the ECL technique.

Performed under reducing conditions.

**Predicted band size:** 88 kDa

**Observed band size:** 88 kDa

**Exposure time:** 10 minutes



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Transferrin Receptor 2/TFR2 antibody (ab80194)

IHC image of Transferrin staining in human normal liver formalin fixed paraffin embedded tissue section, performed on a Leica Bond™ system using the standard protocol F. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH6, epitope retrieval solution 1) for 20 mins. The section was then incubated with ab80194, 5µg/ml, for 15 mins at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors