

Anti-TRIM30/Rpt1 antibody ab76972

1 [References](#) [画像数](#) **2**

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-TRIM30/Rpt1 antibody |
| 製品の詳細 | Rabbit polyclonal to TRIM30/Rpt1 |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: WB, IHC-P |
| 種交差性 | 交差種: Mouse |
| 免疫原 | Synthetic peptide. This information is proprietary to Abcam and/or its suppliers. |
| ポジティブ・コントロール | WB: Mouse spleen tissue lysate. IHC-P: Mouse ovary tissue. |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|----------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー | <p>pH: 7.2</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: PBS</p> |
| 精製度 | Immunogen affinity purified |
| 特記事項(精製) | Affinity chromatography purified via peptide column. |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab76972の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| WB | | Use a concentration of 0.5 - 1 µg/ml. Predicted molecular weight: 57 kDa. |
| IHC-P | | Use a concentration of 10 µg/ml. |

ターゲット情報


関連性

Trans acting factor that regulates gene expression of interleukin 2 receptor alpha chain and of human immunodeficiency virus type 1. May affect IL2R alpha and HIV 1 expression through cis acting negative regulatory elements or through competition with proteins that bind to enhancer or activator sequences.

細胞内局在

Nuclear

画像




Anti-TRIM30/Rpt1 antibody (ab76972) at 1 µg/ml + mouse spleen tissue lysate

Predicted band size: 57 kDa

Observed band size: 57 kDa

Additional bands at: 65 kDa, 70 kDa, 80 kDa. We are unsure as to the identity of these extra bands.



ab76972 at 10 µg/ml staining TRIM30/Rpt1 in mouse ovary tissue section by Immunohistochemistry (Formalin/ PFA fixed paraffin embedded tissue sections).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors