

Anti-TRIOBP antibody [2438C1a] ab58638

画像数 1

製品の概要

製品名	Anti-TRIOBP antibody [2438C1a]
製品の詳細	Mouse monoclonal [2438C1a] to TRIOBP
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Recombinant fragment
免疫原	Recombinant fragment (Human)
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: PBS, 1% BSA
精製度	Protein G purified
特記事項 (精製)	This antibody was purified from culture supernatant of hybridoma cultured in a medium containing bovine IgG depleted (approximately 95%) fetal bovine serum.
ポリ/モノ	モノクローナル
クローン名	2438C1a
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab58638の使用に適用されず**

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration. Predicted molecular weight: 261 kDa.

ターゲット情報

関連性 TRIOBP is a protein that interacts with trio, which is involved with neural tissue development and controlling actin cytoskeleton organization, cell motility and cell growth. This trio binding protein also associates with F actin and stabilizes F actin structures.

細胞内局在 Nucleus. Cytoplasm; cytoskeleton. Note: Localized to F-actin in a periodic pattern.

画像



Anti-TRIOBP antibody [2438C1a] (ab58638) + immunizing recombinant protein

Predicted band size: 261 kDa

Observed band size: 35 kDa

The observed band size does not correspond to the predicted molecular weight as the recombinant immunogen was used for the test lane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors