

Anti-Trypanosoma brucei phosphoglycerate kinase antibody ab89651

★★★★★ [2 Abreviews](#) [画像数 1](#)

製品の概要

製品名	Anti-Trypanosoma brucei phosphoglycerate kinase antibody
製品の詳細	Rabbit polyclonal to Trypanosoma brucei phosphoglycerate kinase
由来種	Rabbit
アプリケーション	適用あり: ICC/IF, WB
種交差性	交差種: Trypanosoma brucei
免疫原	Synthetic peptide derived from the C terminal domain of Trypanosoma brucei phosphoglycerate kinase.

特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituent: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab89651の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
ICC/IF		Use at an assay dependent concentration.
WB	★★★★★ (1)	1/200 - 1/2000. Predicted molecular weight: 55 kDa.

ターゲット情報

関連性

Belongs to the phosphoglycerate kinase family.

画像



All lanes : Anti-*Trypanosoma brucei* phosphoglycerate kinase antibody (ab89651) at 1/200 dilution (in TBS + 5% milk for 12 hours at 4°C)

Lane 1 : Whole cell lysate of *Trypanosoma brucei* - procyclic form

Lane 2 : Whole cell lysate of *Trypanosoma brucei* - bloodstream

Lysates/proteins at 20000000 cells per lane.

Secondary

All lanes : A biotin-conjugated Goat anti-rabbit IgG polyclonal at 1/200 dilution

Predicted band size: 55 kDa

Observed band size: 65,75 kDa

Blocking Step: 5% Milk for 1 hour at 25°C

Detection Method: DAB

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors