

Anti-TUSC5 antibody ab64803

画像数 2

製品の概要

製品名	Anti-TUSC5 antibody
製品の詳細	Rabbit polyclonal to TUSC5
由来種	Rabbit
アプリケーション	適用あり: WB, ICC/IF
種交差性	交差種: Human
免疫原	A synthetic peptide derived from an internal region of human TUSC5
ポジティブ・コントロール	ICC/IF: Hela cells WB: Extracts from LOVO cells
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: PBS, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride
精製度	Immunogen affinity purified
特記事項(精製)	The antibody was affinity purified from rabbit antiserum by affinity chromatography using epitope specific immunogen.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee


Abpromise保証は、次のテスト済みアプリケーションにおけるab64803の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 22 kDa (predicted molecular weight: 19 kDa).
ICC/IF		1/500 - 1/1000.

ターゲット情報

画像



Western blot - Anti-TUSC5 antibody (ab64803)

All lanes : Anti-TUSC5 antibody (ab64803) at 1/500 dilution

Lane 1 : Extracts from LOVO cells

Lane 2 : Extracts from LOVO cells plus peptide at 10µg

Lysates/proteins at 5 µg per lane.

Predicted band size: 19 kDa

Observed band size: 22 kDa



Immunocytochemistry/ Immunofluorescence - Anti-TUSC5 antibody (ab64803)

Immunofluorescence analysis of HeLa cells using ab64803 at a 1:500 dilution.

Left image untreated.

Right image treated with peptide.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors