

### Anti-USP37 antibody ab96244

★★★★★ 1 Abreviews [画像数 2](#)

#### 製品の概要

製品名	Anti-USP37 antibody
製品の詳細	Rabbit polyclonal to USP37
由来種	Rabbit
アプリケーション	<b>適用あり:</b> IHC-P
種交差性	<b>交差種:</b> Human <b>交差が予測される動物種:</b> Chimpanzee, Gorilla 
免疫原	Synthetic peptide, corresponding to a region between C terminal amino acids 930-979 of Human USP37 (NP_065986.1).
ポジティブ・コントロール	Human Breast Carcinoma, Non-Small Cell Lung Cancer, Ovarian Carcinoma and Prostate Carcinoma tissues
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.1% BSA, Tris buffered saline
精製度	Immunogen affinity purified
特記事項 (精製)	ab96244 was affinity purified using an epitope specific to USP37 immobilized on solid support.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション


**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab96244の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/100 - 1/500. Epitope exposure with tris-EDTA pH9.0 buffer is recommended. Likely to work with frozen sections.

## ターゲット情報

組織特異性	Expressed in brain and prostate.
配列類似性	Belongs to the peptidase C19 family. Contains 3 UIM (ubiquitin-interacting motif) repeats.

## 画像



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-USP37 antibody (ab96244)

ab96244, at 1/250 dilution, staining USP37 in formalin-fixed, paraffin-embedded Human ovarian carcinoma by Immunohistochemistry using DAB detection. Epitope retrieval was performed.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-USP37 antibody (ab96244)

ab96244, at 1/250 dilution, staining USP37 in formalin-fixed, paraffin-embedded Human prostate carcinoma by Immunohistochemistry using DAB detection. Epitope retrieval was performed.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors