

Anti-USP40 antibody ab121234

★★★★★ 2 Abreviews    画像数 2

製品の概要

製品名	Anti-USP40 antibody
製品の詳細	Rabbit polyclonal to USP40
由来種	Rabbit
アプリケーション	適用あり: ICC/IF, IHC-P
種交差性	交差種: Human
免疫原	Recombinant fragment, corresponding to amino acids 755-900 of Human USP40.
ポジティブ・コントロール	Human tonsil tissue
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee      Abpromise保証は、 次のテスト済みアプリケーションにおけるab121234の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 0.25 - 2 µg/ml.
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. Use Buffer at pH 6 for antigen retrieval.

## ターゲット情報


機能	May be catalytically inactive.
組織特異性	Broadly expressed.
配列類似性	Belongs to the peptidase C19 family.

## 画像



Immunofluorescent staining of Human cell line U-251MG shows positivity in vesicles and focal adhesions. Recommended concentration of ab121234 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-USP40 antibody (ab121234)



ab121234, at 1/75 dilution, staining USP40 in Human tonsil tissue by Immunohistochemistry using formalin-fixed, paraffin-embedded tissue.

Note: shows strong cytoplasmic positivity in reaction center cells and lymphoid cells outside reaction centra.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-USP40 antibody (ab121234)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors