

Anti-Uteroglobin antibody ab40873

★★★★☆ [6 Abreviews](#) [37 References](#) [画像数 2](#)

製品の概要

製品名	Anti-Uteroglobin antibody
製品の詳細	Rabbit polyclonal to Uteroglobin
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Mouse
免疫原	Recombinant rat uteroglobin His fusion protein.

特記事項

This antibody is useful for distinguishing ciliated from non-ciliated cells in the adult mouse lung.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.01% Sodium azide Constituent: Whole serum
精製度	Whole antiserum
一次抗体 備考	This antibody is useful for distinguishing ciliated from non-ciliated cells in the adult mouse lung.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab40873の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P	★★★★★ (3)	1/1600.

ターゲット情報


機能	Binds phosphatidylcholine, phosphatidylinositol, polychlorinated biphenyls (PCB) and weakly progesterone, potent inhibitor of phospholipase A2.
組織特異性	Clara cells (nonciliated cells of the surface epithelium of the pulmonary airways).
配列類似性	Belongs to the secretoglobin family.
細胞内局在	Secreted.

画像



Adult mouse lung was embedded in paraffin and stained with ab40873 (1:1600 dilution).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Uteroglobin antibody (ab40873)



ab40873 staining Uteroglobin in human lung parenchyma tissue by Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections). Tissue was fixed with formaldehyde and a heat mediated antigen retrieval step was performed using 10mM citrate buffer. Samples were then blocked with 1% BSA, 0.1% Triton-X/ PBS for 15 minutes at 18°C followed by incubation with the primary antibody at a 1/500 dilution for 16 hours at 4°C. An HRP-conjugated pig anti-rabbit polyclonal was used as secondary antibody at a 1/100 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Uteroglobin antibody (ab40873)

Image courtesy of an anonymous Abreview.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors