

Anti-V5 tag antibody ab15828

★★★★★ [1 Abreviews](#) [49 References](#) [画像数 2](#)

製品の概要

製品名	Anti-V5 tag antibody
製品の詳細	Rabbit polyclonal to V5 tag
由来種	Rabbit
アプリケーション	適用あり: ChIP, WB
種交差性	交差種: Species independent
免疫原	Synthetic peptide corresponding to Parainfluenza virus 5 V5 tag aa 50-150 conjugated to keyhole limpet haemocyanin. This ab was raised against ab15829 . Database link: P11207 (Peptide available as ab15829)
ポジティブ・コントロール	WB: Recombinant V5 tagged Groucho homolog. ChIP: A stably transfected 293T human cell line harbouring the GAL4 upstream activation sequence was transiently transfected with a V5 or T7-tagged GAL4 DNA Binding Domain construct.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40 Preservative: 0.02% Sodium azide Constituent: PBS
	Batches of this product that have a concentration < 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our scientific support team who will be happy to help.

精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab15828の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ChIP		Use at an assay dependent concentration. ChIP was performed with 25 ug chromatin, 5ug of antibody and 20 ul of Protein A/G beads.
WB	★★★★★ (1)	Use a concentration of 0.5 µg/ml.

ターゲット情報


関連性 The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids (GKPIPNPLLGLDST), although it has also been used with a shorter 9 amino acid sequence (IPNPLLGLD).

画像



ChIP - Anti-V5 tag antibody (ab15828)

A stably transfected 293T human cell line harbouring the GAL4 upstream activation sequence was transiently transfected with a V5 or T7- tagged GAL4 DNA Binding Domain construct. 48 hours post transfection Chromatin was prepared according to the Abcam X-ChIP protocol. The ChIP was performed with 25 ug chromatin, 5ug of antibody and 20 ul of Protein A/G beads. A non-specific antibody was used as the negative control. The immunoprecipitated DNA was quantified by real time PCR (SYBR Green approach).



All lanes : Anti-V5 tag antibody (ab15828) at 0.5 µg/ml

Lane 1 : Recombinant V5 tagged Groucho homolog

Lane 2 : Recombinant V5 tagged Groucho with V5 tag peptide
(**ab15829**) at 1 µg/ml

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat anti-Rabbit IgG at 1/10000 dilution

Performed under reducing conditions.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors