

Anti-Varicella Zoster Virus gB antibody [10G6] ab52550

画像数 1

製品の概要

製品名	Anti-Varicella Zoster Virus gB antibody [10G6]
製品の詳細	Mouse monoclonal [10G6] to Varicella Zoster Virus gB
由来種	Mouse
アプリケーション	適用あり: ICC/IF, ELISA
種交差性	交差種: Varicella zoster virus
免疫原	Tissue, cells or virus corresponding to Varicella Zoster Virus gB.
ポジティブ・コントロール	ICC/IF: VZV slides
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.40 Constituent: PBS
精製度	Affinity purified
特記事項(精製)	Purity >95% by SDS-PAGE.
ポリ/モノ	モノクローナル
クローン名	10G6
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab52550の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/400.
ELISA		1/102400.

ターゲット情報

関連性 The Varicella zoster virus (VZV) is one of the eight herpes viruses known to affect humans (and other vertebrates). Primary VZV infection results in chickenpox (varicella), which may rarely result in complications including VZV encephalitis or pneumonia.

画像



VZV slides stained with ab52550 at 1/400 dilution followed by anti-Mouse-FITC F(ab')₂ with counterstain

Immunocytochemistry/ Immunofluorescence - Anti-Varicella Zoster Virus gB antibody [10G6] (ab52550)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors